

Perform a Check after completing Setup (Smart Setup or Optional Setup).


How to Perform a Check:

1. You should be in your usual sitting position.
2. Press the Check Button.
3. Refer to the table below.



- Perform a Check frequently, at least once a day!
- If a cover was used while setting up the cushion, the same cover must be used during Check.
- Perform a Check when there is a significant change in altitude.

Symbols and Lights for Check



Check Table

NOTE: After you perform a Check, you may choose to perform a Hand Check to become familiar with your inflation range. See the HAND CHECK GUIDE.

NOTE: You may be able to add or remove a small amount of air and still be at your cushion inflation setting. Smart Check saves your cushion inflation setting with some room for adjustment.

Light	Status	Action
	The cushion is within the saved inflation setting.	No further action is required.
Light	Possible Causes	Action
	Your weight may be unevenly distributed.	Perform a weight shift to circulate air in the cushion, then press the Check Button to get a new reading.
	The quick disconnect is not securely attached to the cushion.	Make sure the connection is secure, then press the Check Button to get a new reading. See TROUBLESHOOTING in this manual.
	The cushion needs more air.	Open the valve and add more air to the cushion (6 to 10 pumps). Close the valve. Press the Check Button to get a new reading. Repeat as needed. If you are unable to make adjustments on your own, discontinue use until assistance is available.
	The cushion or the quick disconnect on the cushion has a slow leak.	See TROUBLESHOOTING in the <i>ROHO® Single Compartment Cushion with Sensor Ready™ Technology Operation Manual</i> , the <i>ROHO® Hybrid Elite® Cushion Operation Manual</i> or contact Customer Care.
Light	Cause	Action
	The cushion has too much air.	Open the valve and remove air from the cushion until you see the green light. Close the valve. If you are unable to perform this task, seek assistance.