

Reference Guide		Customer Design-A-ROHO® Grid														
Cells		Front of Cushion														
4	3	1	2	3	4	5	6	7	8	9	10	11	12	13	14	
2	1	Valve	A	B	C	D	E	F	G	H	I	J	K	L	M	N
Deleted	X	1														
Cut-Out	■	2														
Valve &	↙ ↘	3														
Compartment	1 2	4														
VELCRO	□	5														
		6														
		7														
		8														
		9														
		10														
		11														
		12														
		13														
		14														

### Checklist

Step 1: Specify overall dimensions of cushion (in inches):

\_\_\_\_\_

Step 2: Specify cell heights

\_\_\_\_\_

Proper designation utilized for deleted cells & cut out areas?

Step 3: Specify number of valves

\_\_\_\_\_

Step 4: More than one valve

Compartments indicated?

Step 5: If VELCRO required, please complete

Specify : **Hook or Loop**

Placement of velcro: **Cushion Only**

(Please use dotted lines to indicate) **Cover Only (special)**

**Both**

**Standard Placement**

Step 6: Cover required

Yes / No

Specify Base Model	
HP	___
MP	___
LP	___
QS HP	___
QS MP	___
QS LP	___
ENH	___
CS	___
Other	___

Example DESIGN-A-ROHO: This is a 2R66 with 3" / 7.6 cm cells. The cushion has deleted cells in the back leaving the rubber base. The compartments are divided between rows 3 (C) and 4 (D), and the additional valve is placed in the front right corner.

	1	2	3	4	5	6
Valve	A	B	C	D	E	F
1	3	3	3	3	3	3
2	3	3	3	3	3	3
3	3	3	3	3	3	3
4	3	3	3	3	3	3
5	3	3	X	X	3	3
6	3	3	X	X	3	3

**DO NOT SEND PATIENT PROTECTED HEALTH INFORMATION! IT IS NOT NEEDED TO MAKE THE PRODUCT YOU ARE REQUESTING.**  
Please be advised that custom items take approximately 3-4 weeks from date of order to manufacture.

## Customer Design-A-ROHO® Worksheet

Provider / Distributor Name			
ROHO Account #	Circle one: Quote or Order		
Address			
City	ST	ZIP	Country
Phone	Fax		
Email			
Person Ordering			
Mark For (Reference Name)			
P.O# (Required to place an order)			
Special Instructions: _____			

### How to use the DESIGN-A-ROHO Worksheet

Specific clients may want a ROHO cushion which is designed for their unique needs. This worksheet is intended for that purpose.

The attached grid represents a 14 cell width by 14 cell depth cushion, but cushions with larger dimensions may also be available. Contact ROHO Customer Care for information on any limitations that may exist.

#### Follow these steps when ordering a specialty cushion:

Step 1: Determine the overall size of the cushion required (width by depth) e.g. 15" / 38 cm by 15" / 38 cm

Step 2: The circles represent individual cells. Specify the height of each cell by writing in the desired height. The numbers should be either 1, 2, 3, or 4 which represents 1" / 2.5 cm, 2" / 5 cm, 3" / 7.6 cm and 4" / 10.2 cm, respectively. If you want certain cells deleted, indicate by placing the letter "X" in the circles. When cells are deleted, all that remains is the rubber base. If you want this area removed, please shade in the cell completely on your worksheet. Please note that ROHO will not make a cushion which has only one row of cells on the outer edge of the cushion or one cell standing alone. ROHO suggests that cell variances of 2 or more inches be placed on separate valves for better immersion.

Step 3: Specify the number of valves. If you would like the valve location somewhere other than in the front left corner, please specify "valve" in the area you would like it located.

Step 4: If you have specified more than one valve, indicate compartment divisions by placing a solid line between the compartments.

Step 5: If you would like VELCRO strips attached to the base of the cushion & / or cover, specify the location using dashed lines. Also specify whether you want the hook or loop strip attached. Strips near the edge of the cushion can be both glued and sewn, but strips located toward the center of the cushion can only be glued.

Step 6: Is a cover required?

#### Sizing Chart for HIGH, MID, and LOW PROFILE Cells (All other sizes contact ROHO Customer Care)

1 cell	3" / 7.5 cm	6 cells	11 1/2" / 29.5 cm	11 cells	20" / 50.8 cm
2 cells	4 3/4" / 12 cm	7 cells	13 1/8" / 33.5 cm	12 cells	21 5/8" / 54.9 cm
3 cells	6 3/8" / 16 cm	8 cells	14 7/8" / 38 cm	13 cells	23 3/8" / 59.3 cm
4 cells	8" / 20.5 cm	9 cells	16 1/2" / 42.5 cm	14 cells	25" / 63.5 cm
5 cells	9.75" / 25 cm	10 cells	18 1/4" / 46.3 cm		



100 N. Florida Avenue • Belleville, IL 62221-5429 • USA  
 U.S.: 1.800.851.3449 • 1.618.277.9173 • fax 1.888.551.3449  
 Outside the U.S.: 1.618.277.9150 • fax 1.618.277.6518

[roho.com](http://roho.com)  
[cc@roho.com](mailto:cc@roho.com)