

# ROHO<sup>®</sup> MOSAIC<sup>®</sup> Cushion

- *Conforms to an individual's seated shape*
- *Provides skin/soft tissue protection*
- *Adds little weight to mobility system*
- *Affordable and functional*
- *Adjustable and adaptable to change with the individual over time*



---

*For over 40 years, through research and design, ROHO has manufactured products to fit the needs of a variety of individuals. The MOSAIC cushion is suited for those who require a cushion for skin/soft tissue protection, support and comfort. With AIR FLOATATION™ Technology, the MOSAIC<sup>®</sup> cushion is not only comfortable, but the interconnected air cells distribute forces evenly across the cushion allowing it to be effective for skin/soft tissue protection.*

*A cost effective cushion, the MOSAIC is designed to provide a comfortable, functional and stable sitting environment. MOSAIC is the affordable choice!*

## **COMFORTABLE**

Flexible, interconnected air cells can be adjusted to shape to the body's unique contours, providing customized fit.

## **COST EFFECTIVE**

Constructed of polyvinylchloride, the MOSAIC cushion is an affordable solution for skin/soft tissue protection needs in various facilities or in a home environment.

## **LIGHTWEIGHT**

The MOSAIC cushion weighs only 3/4 pound (.4 kg). Simple to carry from one location to another and easy to use!



# ROHO® MOSAIC® Cushion

## APPLICATIONS

### MOSAIC CUSHION IS APPROPRIATE FOR INDIVIDUALS WHO:

- are at risk for skin/soft tissue breakdown
- currently have skin/soft tissue breakdown
- have normal or impaired sensation
- require accommodation of pelvic asymmetry, up to approximately 1 inch/2.5cm
- require enhanced comfort to their seating surface

## SPECIFICATIONS

### INCLUDED:

MOSAIC cushion, two-way stretch cover, hand inflation pump, repair kit, operation manual and product registration card. Optional Heavy Duty is cover available.

### CONSTRUCTION:

Polyvinylchloride (PVC) made without DEHP

### SIZES:

#### Available to fit chairs

16" x 16"/41cm x 41cm

18" x 16"/46cm x 41cm

18" x 18"/46cm x 46cm

### WEIGHT:

Varies by size, approximately .75lbs/.4kg (based on cushion size 18"x16")

### WEIGHT LIMIT:

250lbs / 113 kg

### WARRANTY:

Cushion: 12-Month Limited Warranty

**U.S. MEDICARE CODE:** E2601

## MORE PRODUCTS FROM ROHO



ROHO® PRODIGY® Mattress Overlay



ROHO® Toilet Seat Cushion



ROHO® AGILITY® Minimum Contour Back System

100 North Florida Avenue • Belleville, IL 62221-5429 • USA • roho.com  
U.S.: 1 800 851 3449 • 1 618 277 9173 • fax 1 888 551 3449  
Outside the U.S.: 1 618 277 9150 • fax 1 618 277 6518

Products in this brochure may be covered by one or more U.S. or worldwide patents. Visit roho.com for further details. Decisions about how to code a product in the U.S. are the responsibility of the provider/supplier. Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local ROHO distributor.

© 2015 ROHO, Inc. The following are registered trademarks of ROHO, Inc: AGILITY®, AIR FLOATATION™, MOSAIC®, PRODIGY®, ROHO®, ROHO, Inc. has a policy of continual product improvement and reserves the right to amend this document. The current version of this document is available on our website at roho.com.

LITD02099 Rev. 12/15 Print 12/15

