

## Performance upgrade

The ROHO® AGILITY® Active Series back supports are offered with an affordable performance carbon fiber upgrade option.

Carbon Fiber Shell

## Move at the speed of lightweight

Starting at only 2.1lbs with integrated hardware, the AGILITY Active series won't weigh you down.

### Specifications

#### ROHO AGILITY Carbon Min

Depth	3"
Width	14", 15", 16", 17" & 18"
Length	10", 13" & 16"
Adjustability	15° of forward & rearward angle adjustment
Hardware	Integrated hardware

#### ROHO AGILITY Aluminum Min

Depth	3"
Width	14", 15", 16", 17" & 18"
Length	8" & 10"
Adjustability	15° of forward & rearward angle adjustment
Hardware	Integrated hardware

#### ROHO AGILITY Carbon Mid

Depth	4"
Width	14", 15", 16", 17" & 18"
Length	10", 13", & 16"
Adjustability	15° of forward & rearward angle adjustment
Hardware	Integrated hardware

#### ROHO AGILITY Aluminum Mid

Depth	4"
Width	14", 15", 16", 17" & 18"
Length	8" & 10"
Adjustability	15° of forward & rearward angle adjustment
Hardware	Integrated hardware

The following are trademarks and registered trademarks of Permobil ROHO Seating and Positioning: AGILITY® and ROHO®.

## ROHO® AGILITY® Active Series

Performance that lasts

**Length options**

Aluminum— 8" & 10"  
Carbon— 10", 13" & 16"

**Width options**

14", 15", 16", 17" & 18"


## Truly functional

Offered in lengths starting at 8" these back supports won't get in the way of what's important. Their compact design and lightweight aluminum or carbon shells allow you to maximize your mobility.

Integrated hardware

Aluminum Shell

15° of forward and rearward angle adjustments



## Hidden strengths

The AGILITY Active Series utilizes durable Zytel discs to achieve the adjustments you need and to provide a rock-solid back support for everyday performance.