



AGILITY Active Carbon Mid Contour Back System

It's made to move you; with compression-molded, lightweight carbon fiber construction, the **AGILITY Active Carbon Mid Contour Back System** combines high-strength performance with clean design. Integrated mounting hardware brings brackets closer to the backrest for streamlined appearance and performance. Plus, the Zytel adjustment disks lock the back in place – but allow for up to 15 degrees of frontward and rearward adjustment.

Options:

Foam with Air Insert

- ✓ Flippable
- ✓ Durable
- ✓ Adjustable AIR FLOATATION insert
- ✓ Strong yet lightweight construction

Dimensions

Width: 13.75 in. to 18.75 in. / 34.9 cm to 47.6 cm

Depth: 4 in. / 10 cm

Height: 10 in. to 16 in. / 25.5 cm to 40.5 cm

Weight: 1.75 lbs. / .815 kg

Based on 14 in. x 10 in. / 35.5 cm x 25.5 cm (weight ranges include hardware)

Weight limit

265 lbs. / 120.2 kg

Back must be properly sized to the individual.

Limited warranty:

Back Shell & Hardware:

36 months

Foam/Cover Components:

18 months

Urethane Inserts: 18 months



Our evidence shows the AGILITY Active Carbon Mid Contour cushion is ideal for people who:

- Require a lightweight and durable back system to increase mobility and daily function
- Require postural support
- Need adjustability to enhance their range of motion and mobility
- Need minimal trunk support