



AGILITY Active Carbon Minimum Contour Back System

The innovative carbon fiber shell is compression molded to provide strength without bulk. The shell can be flipped to accommodate each individual's body type, all while maintaining a higher degree of comfort. With just a single clamp on either side, the AGILITY Active Carbon's integrated hardware maintains a clean appearance while remaining smaller, lighter and stronger than most typical fixed back hardware.

Options:

Foam with Air Insert

- ✓ Flippable
- ✓ Adjustable AIR FLOATATION insert
- ✓ Strong yet lightweight construction
- ✓ Durable

Dimensions

Width: 13.75 in. to 18.75 in. /
34.9 cm to 47.6 cm

Depth: 3 in. / 7.5 cm

Height: 10 to 16 in. /
25.5 cm to 40.5 cm

Weight: 1.75 lbs. / .815 kg

Based on 14 in. x 10 in. /
35.5 cm x 25.5 cm (weight
ranges include hardware).

Weight limit

265 lbs. / 120.2 kg

Back must be properly sized to the
individual.

Limited warranty:

Back Shell & Hardware:

36 months

Foam/Cover Components:

18 months

Urethane Inserts: 18 months



**Our evidence shows the
AGILITY Active Carbon
Minimum Contour cushion
is ideal for people who:**

- Require a lightweight and durable back system to increase mobility and daily function
- Require postural support
- Have good trunk control
- Need adjustability to enhance their range of motion and mobility