



AGILITY Mid Contour Back System

The AGILITY Mid Contour Back System delivers 4 in. (10 cm) of contour depth, to provide posterior trunk support and moderate lateral support. As with all AGILITY Back models, the system is easily installed and adjusted. Height, depth and angle adjustments can be made while you are still in your wheelchair.

Options:

Foam with Air Insert

Foam with Air Insert and Lumbar Support



Optional swing-away and fixed lateral trunk supports



Lateral trunk support



Posterior postural support



Stability and comfort

Dimensions

Width: 13 in. to 20 in. /
33 cm to 51 cm

Depth: 4 in. / 10 cm

Height: 10 in. to 20 in. /
25.5 cm to 51 cm

Weight: 2.5 lbs. to 6 lbs. /
1.1 kg to 2.8 kg

Weight limit

300 lbs. / 136 kg
Back must be properly
sized to the individual.

Limited warranty

Back Shell & Hardware:
5 years

Air Insert: 36 months

Cover & Foam: 18 months

Custom Air Insert: 36 months

**Our evidence shows the
AGILITY Mid Contour
Back System is ideal for
people who:**

- Need minimal trunk support and control
- Require postural support
- Have or have potential for back or shoulder pain