



AGILITY® Minimum Contour Back System

The AGILITY Minimum Contour Back System delivers 2.5 in. (6.5 cm) of contour depth to provide basic postural support in a durable, comfortable package.

AGILITY's unique hardware makes installation and adjustments a breeze. Even while you are in your wheelchair, you can adjust your AGILITY Back height up to 4 in. (10 cm) and your depth up to 1 in. (2.5 cm). You can also get up to 15 degrees of anterior or posterior tilt.

Options:

Foam with Lumbar Support

Foam with Air Insert and Lumbar Support

Foam with Air Insert

Foam with Air Insert and Lumbar Support



Basic posterior postural support



Strength and stability



Optional swing-away and fixed lateral trunk supports



Lightweight and streamlined

Dimensions

Width: 13 in. to 20 in. / 33 cm to 51 cm

Depth: Foam: 2.5 in. / 6.5 cm.; Foam/Air: 3 in. / 7.5 cm

Height: 10 in. to 18 in. / 25.5 cm to 45.5 cm

Weight: 2.5 lbs. to 5.25 lbs. / 1.1 kg to 2.3 kg

Weight Limit

300 lbs. / 136 kg

Back must be properly sized to the individual.

Limited warranty:

Back Shell & Hardware:

5 years

Air Insert: 18 months

Cover & Foam: 18 months

Custom Air Insert: 36 months

Our evidence shows the AGILITY Minimum Contour Back System is ideal for people who:

- Have good trunk control
- Need basic postural support
- Have or have the potential for back or shoulder pain