



ROHO® Custom Cushions

We understand there's no such thing as "one size fits all" when it comes to comfort, safety and health. That's why ROHO® cushions are available in fully customizable shapes, sizes and functionality. Together with your clinician or distributor, you can design your own ROHO – and we'll build it.

- ✓ Skin/soft tissue protection
- ✓ Adjustable
- ✓ Customizable
- ✓ Designed to fit any surface

Dimensions Weight limit: Custom
 Width: Custom
 Height: Custom
 Weight: Custom

Limited warranty
 Cushion: 24 months
 Cover: 6 month



You can:

- Vary the number of valves or compartments to provide side-to-side and front-to-back positioning
- Vary cell heights to provide stability, positioning and skin/soft tissue protection exactly where it's needed
- Add or delete rows of cells beyond dimension of standard size models
- Remove cells or rows
- Add ventilation/drainage holes through base of cushion

Custom Cushions are designed and priced upon request. Restrictions may apply. Please contact Customer Care at cc@roho.com or call (800) 851-3449 or (618) 277-9150. The Design-A-ROHO worksheet is also available on permobilus.com