



Note: cover is required for use.

ROHO® HARMONY® Cushion

The air-filled overlay of the ROHO HARMONY Cushion provides skin/soft tissue protection and positioning, while the contoured foam base delivers just the right amount of stability – particularly on hard surfaces. The HARMONY is a lightweight, comfortable, polyvinyl-chloride cushion that is made without DEHP.

- ✓ Combination of air and foam for ideal seating position
- ✓ Provides comfort and stability
- ✓ Adjustable and lightweight
- ✓ Cost-effective hybrid solution

Dimensions

Width: 13.75 in. to 20 in. /
35 cm to 51 cm

Depth: 16 in. to 17.75 in. /
41 cm to 45 cm

Weight: Approx. 2 lbs. / .9 kg

Weight limit

275 lbs. / 125 kg

Limited warranty

Cushion: 18 months