



## ROHO® HIGH PROFILE® Dual Compartment Cushion

The HIGH PROFILE Dual Compartment Cushion has 4.25 in. (10.5 cm) tall, independent air cells, and each compartment can be adjusted independently to increase positioning and stability for either side-to-side (to adjust for pelvic obliquities) or front-to-back (to adjust bucket or tilt angle). It is designed for those with a distinct pelvic asymmetry or posterior pelvic tilt.

- ✓ Pressure injury prevention and treatment
- ✓ Two independent inflation compartments
- ✓ Aggressive side-to-side and front-to-back positioning
- ✓ Accommodation of pelvic asymmetry, up to approx. 3 in. (7.5 cm)

#### Dimensions

Width: 14.75 in. to 20 in. /  
37.5 cm to 50.5 cm  
Depth: 14.75 in. to 18.25 in. /  
37.5 cm to 46.5 cm  
Height: 4.25 in. / 10.5 cm  
Weight: Approx. 3.25 lbs. /  
1.5 kg

#### Weight limit

No weight limit.  
However, the cushion  
must be properly sized to  
the individual.  
Limited warranty  
Cushion: 24 months  
Cover: 6 months



### Our evidence shows the ROHO HIGH PROFILE Dual Compartment Cushion is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Require accommodation or correction of pelvic asymmetry, approx. 3 in. (7.5 cm)
- Need skin protection, maximum immersion and envelopment
- Have normal or impaired sensation
- Have a history or at risk of skin/soft tissue breakdown