



## ROHO® HIGH PROFILE® Single Compartment Sensor Ready™ Cushion

The tall, interconnected 4.25 in. (10.5 cm) individual air cells envelop you to redistribute pressures and decrease friction and shear. These taller air cells provide a greater contact area and maximum immersion and envelopment – which promotes healing.

With the available Smart Check sensor, you can be sure that with the press of a button, your Sensor Ready cushion is inflated to the proper level.

- ✓ Pressure injury prevention and treatment
- ✓ Adjustability
- ✓ Adaptability
- ✓ Comfort
- ✓ Personal inflation monitoring

#### Dimensions

Height: 4.25 in. / 10.5 cm  
Widths: 11.5 in. to 23.5 in. /  
29 cm to 42 cm

#### Weight limit

No weight limit.  
However, the cushion  
must be properly sized to  
the individual.

#### Limited warranty

Cushion: 24 months  
Cover: 6 months



**Our evidence shows the  
HIGH PROFILE Sensor  
Ready Cushion is ideal  
for people who:**

- Have normal or impaired sensation
- Require accommodation of pelvic asymmetry, up to 3 in. (7.5 cm)
- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)