



## ROHO® LOW PROFILE® Single Compartment Cushion

Designed for those with active lifestyles, the 2.25 in. (5.5 cm) tall, interconnected air cells maximize your mobility for transfers. The cushion provides the right balance of immersion and stability, protecting you from tissue breakdown by redistributing your weight across the cushion.

The LOW PROFILE cushion is also ideal for use as a secondary cushion on an office chair, in the car, or in the bath – anywhere you need extra cushioning and protection from injury.

- ✓ **Lightweight**
- ✓ **Transfer stability**
- ✓ **Adjustable**

#### Dimensions

Width: 11.5 in. to 23.5 in. /  
29 cm to 59.5 cm

Depth: 11.5 in. to 18.25 in. /  
29 cm to 46.5 cm

Height: 2.25 in. / 5.5 cm

Weight: Approx. 2.5 lbs. /  
1.1 kg

#### Weight limit

No weight limit.  
However, the cushion  
must be properly sized to  
the individual.

#### Limited warranty

Cushion: 24 months  
Cover: 6 months



### Our evidence shows the LOW PROFILE Single Compartment Cushion is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Have normal or impaired sensation
- Have a history or are at risk of skin/soft tissue breakdown