



ROHO® MID PROFILE™ Single Compartment Cushion

The 3.25 in. (8.5 cm) individual air cells bring you the perfect balance immersion, protection and positioning control, redistributing pressure to decrease friction and shear. With a simple inflation adjustment, your ROHO conforms to your shape, minimizing tissue deformation and maximizing blood flow.



Pressure injury prevention and treatment



Lightweight



Adjustable



Washable



Stability

Dimensions

Width: 13.25 in. to 20 in. /
33.5 cm to 50.5 cm

Depth: 13.25 in. to 20 in. /
33.5 cm to 50.5 cm

Height: 3.25 in. / 8.5 cm

Weight: Approx. 3.0 lbs. /
1.4 kg

Weight limit

No weight limit.
However, the cushion
must be properly sized to
the individual.

Limited warranty

Cushion: 24 months
Cover: 6 months



Our evidence shows the MID PROFILE Single Compartment Cushion is ideal for people who:

- Currently have, have a history of, or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Need skin protection, immersion and envelopment
- Have normal or impaired sensation