



ROHO® MID PROFILE™ Single Compartment Sensor Ready™ Cushion

The 3.25 in. (8.5 cm) individual air cells provide the ideal level of immersion and skin protection, thanks to the Shape Fitting Technology® of ROHO DRY FLOATATION®. The shape and depth of the cells adapt to your movements and conform to your body's changes throughout the day.

With the addition of the Smart Check device, you gain additional peace of mind, knowing your ROHO is always inflated to the proper level.



Pressure injury prevention and treatment



Adjustability



Adaptability



Comfort



Personal inflation monitoring

Dimensions

Widths: 13.25 in. to 20 in. /

34 cm to 51 cm

Height: 3.25 in. / 8 cm

Weight: Approx. 3 lbs. / 1.4 kg

Weight limit

No weight limit.

However, the cushion

must be properly sized to the individual.

Limited warranty

Cushion: 24 months

Cover: 6 months



Our evidence shows the MID PROFILE Sensor Ready Cushion is ideal for people who:

- Have normal or impaired sensation
- Require accommodation of pelvic asymmetry, up to approximately 2 in. (5 cm)
- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)