

ROHO® MINI-MAX® Cushion

This is a cushion designed for short-term use during recreational and sport activities, bathing, traveling, etc. The MINI-MAX provides shock and impact absorption, as well as skin protection. The 1.25 in. (3.2 cm) air-filled cells are fully adjustable.



- ✓ Comfortable
- ✓ Absorbs shock
- ✓ Adjustable

Dimensions

Height: 1.25 in. / 3 cm
Width: 15 in. to 18.25 in. /
38.5 cm to 46.5 cm
Depth: 14.75 in. to 16 in. /
37.5 cm to 40.5 cm
Weight: 2 lbs. / 0.9 kg

Weight limit

None, if properly sized

Limited warranty

Cushion: 24 months
Cover: 6 months