

ROHO® PACK-IT® Cushion

The ROHO PACK-IT Cushion is a smaller version of the MINI-MAX. It's designed for athletic competitions, bathing, traveling or as a scapular, sacral or lumbar support. The PACK-IT's size and built-in strap makes it easy to attach and remove the cushion on the go.



Strap for easy attachment and removal



Absorbs shock



Skin/soft tissue protection



Adjustable



Dimensions

Height: 1.25 in. / 3 cm
Width: 16.25 in. / 41.5 cm
Depth: 9.5 in. / 24 cm
Weight: 1.5 lbs. / 0.7 kg

Weight limit

None, if properly sized.
Limited warranty
Cushion: 24 months
Cover: 6 months