



## ROHO PRODIGY Mattress Overlay System

The PRODIGY's three independent sections can be adjusted separately to accommodate your specific, changing needs. The flexible, interconnected AIR FLOATATION air cells match your body's contour to provide skin/soft tissue and deep tissue injury protection.

The modular mattress sections can be adjusted to a flat base or on that is anatomically shaped. The PRODIGY Mattress Overlay System is not dependent on an electrical air pump and thus it produces no sound, heat or vibrations.

The PRODIGY Mattress Overlay has drainage holes in the base, allowing fluid to drain away from the body.

- ✓ Adjustable to size and weight change over time
- ✓ Assists with pain management
- ✓ Provides 3 in. (7.5 cm) of immersion
- ✓ Easy to clean, disinfect and patch
- ✓ Polyvinyl chloride (made without DEHP)

#### Dimensions

Height: 3.5 in. / 9 cm

Width: 36 in. / 91.5 cm

Depth: 81.75 in. / 207.5 cm

Weight: 9.5 lbs. / 4.3 kg

#### Weight limit

300 lbs. / 136 kg

The overlay must be properly sized to the individual and bed mattress.

#### Limited warranty

12 months