



# ROHO® QUADTRO SELECT® HIGH PROFILE® Cushion

The QUADTRO SELECT® HIGH PROFILE® features four distinct compartments, each controlled by our exclusive ISOFLO® Memory Control. This gives you ultimate control over the inflation levels across your ROHO cushion, making it easy to ensure you are positioned properly. The interconnected 4.25 in. (10.5 cm) tall air cells envelop your body, delivering maximum immersion and positioning, while protecting you from skin and soft tissue breakdown.



**Pressure injury prevention and treatment**



**Stability**



**In-place customizability with ISOFLO® Memory Control**



**Maximum positioning control**



**Four independent compartments for individualized fit**

#### Dimensions

Width: 12.75 in. to 26.25 in. / 32.5 cm to 67 cm

Depth: 11.75 in. to 67 in. / 30 cm to 51.5 cm

Height: 4.25 in. / 10.5 cm

Weight: Approx. 3.5 lbs. / 1.6 kg

#### Weight limit

No weight limit. However, the cushion must be properly sized to the individual.

#### Limited warranty

Cushion: 36 months  
Cover: 6 months



## Our evidence shows the QUADTRO SELECT HIGH PROFILE is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Require accommodation or correction of pelvic asymmetry, up to approx. 3 in. (7.5 cm)
- Require pelvic and lower extremity positioning
- Have asymmetrical posture and need front-to-back and side-to-side positioning
- Need additional stability while seated
- Tend to lean or slide in their wheelchair