



ROHO® QUADTRO SELECT® MID PROFILE™ Cushion

The QUADTRO SELECT® MID PROFILE™ features four distinct compartments, each controlled by our exclusive ISOFLO® Memory Control. This gives you ultimate control over the inflation levels across your ROHO cushion, making it easy to ensure you are positioned properly. The interconnected 3.25 in. (8.5 cm) air cells deliver the ideal combination of customizable positioning and stability, all while protecting you from skin and soft tissue breakdown.

- ✓ Pressure injury prevention and treatment
- ✓ Four independent compartments for individualized fit
- ✓ In-place individualization with ISOFLO® Memory Control
- ✓ Blend of positioning and stability

Dimensions

Width: 14.5 in. to 19.5 in. /
36.5 cm to 49.5 cm

Depth: 13.5 in. to 18.5 in. /
34.5 cm to 47 cm

Height: 3.25 in. / 8.5 cm

Weight: Approx. 3 lbs. /
1.3 kg

Weight limit

No weight limit.
However, the cushion
must be properly sized to
the individual.

Limited warranty

Cushion: 36 months
Cover: 6 months



Our evidence shows the QUADTRO SELECT MID PROFILE is ideal for people who:

- Currently have, have a history of, or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Require pelvic and lower extremity positioning
- Have asymmetrical posture and need front-to-back and side-to-side positioning