

ROHO® Shower/Commode Cushion

Using ROHO's exclusive DRY FLOATATION® Technology, the adjustable, air-filled Shower/Commode cushion provides much-needed cushioning and protection for bedside commodes and shower wheelchairs. Designed for use without a cover, the cushion is easy to clean and disinfect.

- ✓ Adjustable
- ✓ Comfortable
- ✓ Easy to clean and disinfect
- ✓ Skin/soft tissue protection



Dimensions

Height: 2.25 in. / 5.5 cm
Width: 14.75 in. / 37.5 cm
Depth: 14.75 in. / 37.5 cm
Weight: 1.5 lbs. / 0.7 kg

Weight limit

None, if properly sized
Limited warranty
24 months