

ROHO® Toilet Seat Cushion

The adjustable, air-filled cellular neoprene surface is designed for short-term use on toilets and portable bedside commodes. It snaps easily into place and is easy to clean and disinfect.

- ✓ Quickly snaps into place
- ✓ Adjustable
- ✓ Easy to clean and disinfect
- ✓ Skin/soft tissue protection



Dimensions

Height: 2.25 in. / 5.5 cm
Width: 7.75 in. / 20 cm
Depth: 15.75 in. / 40 cm
Weight: 1.25 lbs. / 0.5 kg

Weight limit

None, if properly sized
Limited warranty
24 months