



Note: cover is required for use.

## ROHO® nexus SPIRIT® Cushion

The nexus SPIRIT Cushion combines ROHO DRY FLOATATION® Technology with the stability of a soft, contoured foam base. This base provides stability and positioning and support when transferring. The individual air cells envelop your body, providing pressure relief and preventing skin/soft tissue breakdown.

- ✓ Combination of air and foam for positioning
- ✓ Soft contoured foam base for stability
- ✓ Transfer stability and support
- ✓ Pressure injury prevention and treatment
- ✓ Accommodates up to approx. 1 in. (2.5 cm) pelvic asymmetry

### Dimensions

Width: 14 in. to 20 in. /  
39.5 cm to 49.5 cm

Depth: 14 in. to 18 in. /  
38 cm to 36 cm

Air Insert Height: 2.25 in. /  
5.5 cm

Total Height: 3.5 in. / 9 cm

Weight: Approx.  
3 lbs. / 1.3 kg

### Weight limit

No weight limit.  
However, the cushion  
must be properly sized to  
the individual.

### Limited warranty

Cushion: 24 months  
Cover: 6 months