



Superstand HLT

Pediatric Standing System

Canadian Order Form/Price List, effective January 2020



Complete Positioning



Easy Loading



Portability

All Prime products are made to order and are non-returnable and non-exchangeable.



Client: _____ Age: _____
 Height: _____ Weight: _____ Dx: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____
 Facility: _____
 Therapist: _____ Phone: _____

PO: _____ Date: _____
 Dealer: _____
 Contact: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____ Fax: _____
 E-mail: _____

Standard Features: Horizontal high-loading surface; Prone, Vertical, Supine or Multi-Position configurable; hip and chest pads; hands-free standing angle adjustments; height and depth adj knee system; lock out block and 3" locking casters.

Specifications subject to change without notice.

Superstand HLT EI
 Client Height: Supine 20" - 38" Prone 22" - 40"
 Base: 19" wide x 28" long. Loading Height 27"
 Client weight limit 65 pounds

Superstand HLT Size 1
 Client Height: Supine 28" - 48" Prone 30" - 50"
 Base: 21" wide x 37 3/4" long. Loading Height 28"
 Client weight limit 125 pounds

Superstand HLT Size 2
 Client Height: Supine 40" - 60" Prone 40" - 62"
 Base: 25" wide x 45" long. Loading Height 29"
 Client weight limit 150 pounds

Folding Base: Includes quick release folding base, mast mounting hardware, gas cylinder and 4 plastic locking casters.
 Base Size: **EI \$1,457.50** **Size 1 \$1,457.50** **Size 2 \$1,787.50**
 Unit Configuration / Setup: **Supine** **Prone** **Multi-Position** **Custom**
 Color: **Candy Red** **Candy Blue** **Candy Purple** **Neon Orange** **Neon Pink** **Shocker Green** **Black** **White** **Razzberry**

Mast
 EI Mast: EI Mast Assembly fits 20" - 38" client, includes chest and hip support pad knee supports and solid foot platform.....\$ **2,161.50**
 Size One Mast: Assembly fits 28" - 48" client, includes chest and hip support pad knee supports and solid foot platform.....\$ **2,161.50**
 Size Two Mast: Assembly fits 40" - 60" client, includes chest and hip support pad knee supports and solid foot platform.....\$ **2,662.00**
 ***Note: EI Mast fits in all bases, Size One mast fits in Size one and Size Two base, Size Two mast fits in Size Two base only.
 Upholstery Color: **Black Vinyl** **Blue Vinyl** **Black** **Red Vinyl** **Custom**

Lateral Positioning Supports

<input type="checkbox"/> Lateral Chest Supports EI: 3 1/2" deep by 3" high.....	\$	253.00
<input type="checkbox"/> Lateral Chest Supports Size 1: 4 1/2" deep by 4" high.....	\$	253.00
<input type="checkbox"/> Lateral Chest Supports Size 2: 4 1/2" deep by 6" high.....	\$	291.50
<input type="checkbox"/> Lateral Hip Supports EI: 3 1/2" deep by 3" high.....	\$	253.00
<input type="checkbox"/> Lateral Hip Supports Size 1: 4 1/2" deep by 5 1/2" high.....	\$	253.00
<input type="checkbox"/> Lateral Hip Supports Size 2: 4 1/2" deep by 6 1/2" high.....	\$	291.50
<input type="checkbox"/> Sliding Adj. Tracks on Size 1 Non-Contoured Chest pad or on Size 2 Standard Chest pad***:	\$	239.80
<input type="checkbox"/> Sliding Adj. Tracks on Hip pad***: Sliding tracks supply maximum adj. for lateral pads for height and angle.....	\$	239.80

***Note: Sliding adjustment tracks not available on EI Unit.

Chest Pad lateral mounting brackets: 90° Standard 1" Offset 2" Offset
Hip Pad lateral mounting brackets: 90° Standard 1" Offset 2" Offset

EI Lateral Support Width Adjustment range 90° 6" - 10" ~ 1" offset 4" - 8" ~ 2" offset 3" - 6"	Size 1 Lateral Support Width Adjustment range 90° 8" - 12" ~ 1" offset 6" - 10" ~ 2" offset 4" - 8"	Size 2 Lateral Support Width Adjustment range 90° 10" - 14" ~ 1" offset 9" - 12" ~ 2" offset 6" - 11"
--	---	---

Straps and Supports

<input type="checkbox"/> Chest Belt: <input type="checkbox"/> EI <input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 Black nylon velcro chest belt.....	\$	94.60
<input type="checkbox"/> Padded Chest Belt: <input type="checkbox"/> EI <input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 Black nylon chest belt with pad.....	\$	171.60
<input type="checkbox"/> Body Support Vest: <input type="checkbox"/> EI <input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 Quick release vest for upper body support.....	\$	288.20
<input type="checkbox"/> Soft Hip Support: <input type="checkbox"/> EI <input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 Soft strap for hip support while in Prone.....	\$	245.30

Harness Colors: **Black Denim** **Blue Denim** **Black Vinyl** **Blue Vinyl** **Red Vinyl** **Customer supplied** **Custom**

Head Support:

<input type="checkbox"/> Headrest: Height and depth adjustable padded headrest (includes 1/2" Square Tube Mount)	\$	385.00
<input type="checkbox"/> 1/2" Square Tube Mount: Height and depth adjustable mount for any headrests that use a 1/2" horizontal mounting bar.....	\$	171.60
<input type="checkbox"/> i2i Flip Back Mount: Flip back mounting hardware for the i2i headrest (hardware only)	\$	462.00
<input type="checkbox"/> 1/2" Square Tube Flip Back Mount: Flip back mount for any headrest using a 1/2" horizontal mounting bar.....	\$	319.00
<input type="checkbox"/> Chinrest: Height adjustable chin support.....	\$	171.60
<input type="checkbox"/> Padded Head Strap: Sheepskin-lined strap with Velcro closure.....	\$	165.00

All Prime products are made to order and are non-returnable and non-exchangeable.





Lower Extremity Positioning

- Independent Footplates:** Individually height adjustable footplates with heel loops.....\$ **332.20**
- Multi-Position Foot System:** Individually adj. for planter-flexion, dorsi-flexion, rotation, height, with heel loops and foot straps...\$ **621.50**
- Size-Adjustable Knee Pads (pair):** Allows knee pad width adjustment from 2 1/2" - 4 1/2" in width (not available in EI).....\$ **478.50**
- Contracture Knee System:** Supports contracted Lower Extremities with a fully adjustable system (on System)\$ **1,060.97**
- D-Ring Footstraps:** 1" Velcro D-Ring straps.....\$ **94.60**
- Shoe Holders (pair):** Formed ABS shoe holders with ankle and toe straps.....\$ **275.00**
 Size: Sm (6" long x 4 1/4" wide) Med (8" long x 4 7/8" wide) Lrg (10 1/2" long x 6" wide)
- Height Adjustable Pommel:** Assists in maintaining lower extremity position.....\$ **319.00**
- Pommel W/ Rigid Hip Pad:** Height adjustable pommel with height and depth adjustable rigid hip pad.....\$ **572.00**
- ABduction Leg System:** Includes rotating knees and fully adjustable shoe holders.....\$ **649.00**
- ABduction Leg System with Pommel:** Includes rotating knees, fully adj shoe holders and pommel with rigid pad.....\$ **866.80**

Trays

- Opaque Upper Extremity Positioning Tray:** Multi-Positional tray with height, depth and angle adjustment.....\$ **691.90**
- Clear Upper Extremity Positioning Tray:** Multi-Positional tray with height, depth and angle adjustment.....\$ **861.30**
- Tray Options:** Opaque overlay.....\$ **231.00** Clear overlay.....\$ **231.00** Tray Easel.....\$ **345.50** Edging.....\$ **138.60**
- Activity Bowl:** Tray cut-out with 7" bowl (clear tray only)..... Add \$ **239.80**
- Padded Tray Covers:** Full size removable padded cover.....\$ **283.80** Wrap around non removable elbow padding.....\$ **283.80**

Casters

- Steel Caster Upgrade:** High performance easy rolling steel casters that lock for both swivel and roll.....\$ **275.00**



Standard Superstand HLT Support Pad Specifications:

Specifications subject to change without notice.

EI: Chest Pad: 8 1/2" W x 7" H ~ Chest Laterals: 3 1/2" D x 3" H ~ Hip Pad: 8 1/2" W x 5" H ~ Hip Laterals: 3 1/2" D x 3" H ~ Knee Pads: 2" W x 3" H
Size 1: Chest Pad: 10" W x 8 1/2" H ~ Chest Laterals: 4 1/2" D x 4" H ~ Hip Pad: 10" W x 8" H ~ Hip Laterals: 4 1/2" D x 5 1/2" H ~ Knee Pads: 4" W x 3 1/2" H
Size 2: Chest Pad: 12 1/2" W x 15" H ~ Chest Laterals: 4 1/2" D x 6" H ~ Hip Pad: 12 1/2" W x 9" H ~ Hip Laterals: 4 1/2" D x 6 1/2" H ~ Knee Pads: 4 1/4" W x 5" H
Headrest Specifications: **EI** and **Size 1** 6 1/2" W x 3 1/2" H **Size 2** 8" W x 4 1/2" H

Custom Configurations: We are proud of our proven commitment to serving the needs of our customers with the finest standing products available since 1984. Please call us when you need something special for your client. Supply the information in the following area so we can know your needs and we will assist in every way that we can.


Client Measurements:

A: Thorax Width ___" B: Hip Width ___" C: Hip-to-Knee ___" D: Knee-to-Heel ___" E: Shoe Size ___" F: Shoulder Height ___" G: Height ___"
 H: Inseam ___" I: Weight ___lbs. **Note:** All Measurements Assume Symmetrical Standing and/or Sitting Posture.

Custom Pads: Use this area for custom-sized pads on your Superstand HLT. All custom pads must be quoted for price and manufacturing time.

Chest Pad: ___" Wide x ___" High Standard Shape Pad with Top Contour Custom Shape, describe below ~ Chest Laterals ___" Wide x ___" High
 Hip Pad: ___" Wide x ___" High Standard Shape Pad Custom Shape, describe below ~ Hip Laterals ___" Wide x ___" High

Custom Tray Surface: ___" Wide X ___" Deep; Standard Tray is 22 1/2"w x 19"d for Size one, and 28"w x 24"d for Size two.



Include information for any tray cutouts or modifications desired in this area or submit a paper or cardboard template. Other information may be added in the area below.

Special Requests / Additional Information: _____

All Prime products are made to order and are non-returnable and non-exchangeable.

