



Build to Suit Multi-Position Pediatric Standing Frame

BTS Multi-Position Stander & BTS Youth Multi-Position Stander

Canadian Order Form/Price List, effective January 2020

Client: _____ Age: _____
 Height: _____ Weight: _____ Dx: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____
 Facility: _____
 Therapist: _____ Phone: _____

PO: _____ Date: _____
 Dealer: _____
 Contact: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____ Fax: _____
 E-mail: _____

Multi-Position Build To Suit Pediatric Standing Frame

Standard Features: Prone, Vertical and Supine capable; Hip support pad with lateral supports; Upper body support pad; Velcro chest belt; Manual standing angle adjustment from 30 to 90 degrees in 7 degree increments; Height, depth and width adjustable knee system; One-piece depth adjustable foot platform with heel block and 3" plastic casters.

<input type="checkbox"/> BTS Multi-Position Stander: Client Ht. 30" - 48" With Standard Foot System 28" - 48" with Height Adjustable or Multi- Position Foot Systems Base: 21" wide x 40" deep \$2,139.50		<input type="checkbox"/> BTS YOUTH Multi-Position Stander: Client Ht. 43" - 60" With Standard Foot System 40" - 60" with Height Adjustable or Multi- Position Foot System Base: 24" wide x 48" deep \$3,140.50
--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Frame Color: <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Candy Purple Support Vest: <input type="checkbox"/> Black Denim <input type="checkbox"/> Blue Denim	Upholstery Color: <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Black <input type="checkbox"/> Black Vinyl <input type="checkbox"/> Blue Vinyl <input type="checkbox"/> Red Vinyl
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Lateral Brackets: **BTS** is standard with 5 1/2"H X 4"D Hip Laterals with 90 Degree Brackets (width adjustment from 9" - 13").
BTS Youth is standard with 6 1/2"H X 4 1/2"D Hip Laterals with 90 Degree Brackets (width adjustment from 11" - 15").
 1" Offset Upgrade \$33.00 Width adjustment of 7" - 11" on BTS, and 9" - 13" on BTS Youth
 2" Offset Upgrade \$33.00 Width adjustment of 5" - 9" on BTS, and 7" - 11" on BTS Youth

Hip Brackets: 90 Degree 1" Offset 2" Offset **Chest Brackets:** 90 Degree 1" Offset 2" Offset

Body Support Options:

Lateral Chest Supports: Pair of lateral chest pads.....\$ 203.50
 BTS MPS: 3-1/2" x 3-1/2" pads with 90 degree brackets.
 BTS MPS Youth: 6-1/2" x 4-1/2" pads with 90 degree brackets.

Upper Body Support Vest: Quick release vest for upper body positioning.....\$ 185.90

Padded Head Strap:\$ 115.50

Headrest: Ht and depth adj. padded headrest.....\$ 225.50

Chinrest: Height adjustable chin support.....\$ 141.90

Frame Options:

Steel Caster Upgrade: Heavy duty casters lockout swivel and roll.....\$ 236.50

Pneumatic Assisted Standing Angle Adjustment:.....\$ 247.50

Lower Extremity Positioning Options:

Pommel: Ht-adj. pommel assembly..... \$ 247.50

Pommel w/ Ht. Adj. Rigid Hip Pad:\$ 533.50

Size Adjustable Knee Pads: Allow for individual knee pad width adjustment of 2-1/2" - 4-1/2" for growth.....\$ 280.50

Foot Options:

Heel Loops: 2" high web heel loops..... \$ 60.50

Footstraps: 1" Velcro straps..... \$ 60.50

D-Ring Footstraps: 1" Velcro D-Ring straps..... \$ 93.50

Independent Footplates: Individually ht. adj. footrests\$ 181.50

Multi-Position Foot System: Individually adjustable for planter-flexion, dorsi-flexion, rotation and height\$ 418.00

Shoe Holders (pair): ABS shoe holders with ankle and toe straps

Sm (6" long x 4-1/4" wide) Med (8" long x 4-7/8" wide)

Lrg (10-1/2" long x 6" wide)\$ 236.50

Tray Options:

Tray: Height, depth and angle adjustable opaque tray.....\$ 555.50

Clear Tray Upgrade: Clear tabletop upgrade..... \$ 137.50

Tray Overlay:\$ 198.00

Clear overlay Opaque overlay

Activity Bowl: Tray cut-out with 7" bowlAdd \$ 176.00

Padded Tray Cover: Removable padded cover\$ 247.50

Tray Easel: Free standing easel\$ 345.40

All Prime products are made to order and are non-returnable and non-exchangeable.



Distributed in Canada by Permobil www.permobil.ca
 Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944


Build to Suit Multi-Position Pediatric Standing Frame

BTS Multi-Position Stander & BTS Youth Multi-Position Stander

Canadian Order Form/Price List, effective January 2020

Client Name: _____

Page 2



Custom Support Pads: Use this area if you need custom support pads on your **BTS**. Custom pads must be quoted for price and manufacturing time.

Hip Pad: ____" Wide x ____" High Standard size is 10" wide x 7" high for **BTS** and 12" wide x 8" high for **BTS Youth**.


Hip Laterals ____" Wide x ____" High Standard size is 4" deep x 5 1/2" high for **BTS** and 4 1/2" deep x 6 1/2" high for **BTS Youth**.

Chest Pad: ____" Wide x ____" High Standard size is 10" wide x 8" high for **BTS** and 12" wide x 15" high for **BTS Youth**.


Chest Laterals ____" Wide x ____" High Standard size is 3 1/2" deep x 3 1/2" high for **BTS** and 4 1/2" deep x 6 1/2" high for **BTS Youth**.

Custom Tray Surface:

Tray: ____" Wide x ____" Deep Standard Tray is 22 1/2"w x 19"d for **BTS** and 28"w x 24"d for **BTS Youth**.



Please supply any information for cutouts or modifications desired in this area or submit a paper or cardboard template. Additional information may be added in the special needs area below.



Client Measurements:

A: Thorax Width _____" B: Hip Width _____" C: Hip-to-Knee _____"
 D: Knee-to-Heel _____" E: Shoe Size _____" F: Shoulder Height _____"

Note: All Measurements Assume Symmetrical Standing and Sitting Posture.

Special Needs:

We are proud of our proven commitment to serving the needs of our customers with the finest standing products available since 1984. We encourage you to call us with your special needs, questions, or any item that we may assist with in your standing program. Our standing and positioning specialists will be happy to answer your questions or assist you with the design and fabrication of a unit that meets your needs. Please supply as much of the information in the client measurement area of this form as possible, and list special needs in the following area to speed up the design process and / or response time.

All Prime products are made to order and are non-returnable and non-exchangeable.



Distributed in Canada by Permobil www.permobil.ca
 Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944