



Granstand III & Kidstand III

Modular Hydraulic Standing System

Canadian Order Form/Price List, effective January 2020

Client: _____ Age: _____

Height: _____ Weight: _____ Dx: _____

Address: _____

City, Prov, PC: _____

Phone: _____

Facility: _____

Therapist: _____ Phone: _____

PO: _____ Date: _____

Dealer: _____

Contact: _____

Address: _____

City, Prov, PC: _____

Phone: _____ Fax: _____

E-mail: _____



Granstand III Client Height: 5' 6" - 6' 5"

\$4,218.00 Weight Limit 250 lbs

Granstand III Client Height: 4' 6" - 6' 5"
Requires Mobile Foot System

Granstand III / Kidstand III Modular Standing Systems include our **NO TRANSFER DESIGN**, allowing the client to stand directly from their mobility product. Granstand III base design lifts from the wheelchair, bed, mat table or other surfaces and re-seats client safely and properly onto any of these surfaces.

Kidstand III Client Height: 3' 6" - 4' 8"

\$3,712.50 Weight Limit 250 lbs

Kidstand III Client Height: 3' - 4' 8"
Requires Mobile Foot System



Standard Features: Expandable transport base with 3" plastic casters; Height and depth adjustable chest pad; Height and depth adjustable clear tabletop; Width and depth adjustable foot pads; **Pull-to-Place** seat sling; and Depth and height adjustable knee pad.

Pump Handle: Right Left Frame Colors: Black Candy Red Candy Blue Candy Purple Shocker Green
 Neon Orange Neon Pink White Razzberry

Granstand III Mobile Kit
Requires Mobile Foot System

\$ 1,281.50

Kidstand III Mobile Kit
Requires Mobile Foot System

The **Mobile Kit** allows the client with upper extremity function to be mobile while standing. The Granstand III / Kidstand III MSS with a mobile package allows it to be used as a self-propelled lift / transfer device, as well as a complete standing system.

Each mobile system includes handrim locks, handrims, steel caster upgrade and front drive wheels.



Seat Sling Options:

- Sheepskin Sling Cover:** Genuine sheepskin cover with Velcro attachments for standard web sling.....\$137.50
 - Pocket Seat Sling Cover:** Pocket slipcover design to allow for additional padding of the standard web sling..... \$209.00
 - Vinyl Seat Sling:** One piece vinyl **pull-to-place** seat sling allows greater surface contact and ease of cleaning.....\$236.50
 - Contoured Vinyl Sheepskin Lined Seat Sling:** One piece vinyl **pull-to-place** seat sling with full sheepskin lining and mild contour.....\$280.50
- Sling Sizes: Kidstand Sling 7" H x 19" W (std) Granstand Slings Med 9" H x 22"W Large 11"H x 24"W (std), XL 13.5"H x 29"W

If no sling size is marked the standard sling for the product will be included.

Lower Extremity Positioning Options:

- Standard Foot Pads:** Width and depth adjustable. (Products are not mobile in use with this item).....No Charge
- Height Adjustable Mobile Foot System:** Adjusts in width, depth, height, toe in and toe out. Matches the foot positioners to the height of the wheelchair footplates. Supplies maximum growth potential, and allows the Granstand III to be mobile while in use.....\$533.50
- Multi Adjustable Mobile Foot System:** Adjusts in width, depth, height, toe in, toe out and individual planter-flexion or dorsi-flexion positioning. Supplies maximum growth potential and height adjustment, plus allows the Granstand III to be mobile while in use.....\$533.50
- Multi-Adjustable Knee System:** Tool-free mechanism easily adjusts in height, width and depth for optimum lower body positioning. Independent width positioning allows for proper placement of the knee pad and fast adjustment for multiple client use.....\$522.50
- ROHO Knee Pads:** Additional padding for clients with special positioning needs such as contractures (attaches to either knee system).....\$500.50
- Foot Straps:** Velcro straps for assistance in foot positioning.....\$64.90
- D-Ring Foot Straps:** Velcro straps with D-Rings for additional holding strength.....\$93.50
- Foot Risers:** Available in 1.5" high increments (*for standard foot pads only*).....\$181.50



Standard Foot Pads



Ht. Adj. Foot System



Multi-Adj. Foot System



Std. Knee Pad



Multi-Adj. Knee System



Roho Knee Pads

All Prime products are made to order and are non-returnable and non-exchangeable.



Distributed in Canada by Permobil www.permobil.ca
Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944


Granstand III & Kidstand III

Modular Hydraulic Standing System
Canadian Order Form/Price List, effective January 2020

Client Name: _____

Upper Body Positioning Options:

- Padded Chest Belt.....\$ 181.50
- Small Body Pad (4" x 6").....\$ 55.00
- Contoured Body Pad.....\$ 236.50
- Contoured Body Pad with Lateral Supports.....\$ 495.00
- Contoured Body Pad with Laterals and Padded Support Belt...\$ 566.50
- Custom Body Pad.....Call for Quote



Base / Frame Options:

- Steel Caster Upgrade.....\$ 247.50
- Quad Pump Handle.....\$ 181.50
- 3" Wider Frame.....\$ 852.50
- Two Piece Pump Handle.....\$ 137.50




Tray Options:

- Pneumatic Tray Lift.....\$ 319.00
- Tray Easel.....\$ 324.50
- Tray Edging.....\$ 134.50


Tray Overlays and Covers:

- Opaque Tray Overlay.....\$ 181.50
- Clear Tray Overlay.....\$ 181.50
- Padded Tray Cover.....\$ 286.00

Custom Tray Surface: ___" Wide x ___" Deep (Standard Tray is 24"w x 18"d)



Please supply any information for cutouts or modifications desired in this area or submit a paper or cardboard template. Additional information may be added in the special needs area below.



Call for Quote

Client Measurements:

A: Thorax Width _____" B: Hip Width _____" C: Hip to Knee _____"
 D: Knee to Heel _____" E: Shoe Size _____" F: Shoulder Height: _____"

Note: All Measurements Assume Symmetrical Standing and Sitting Posture.

Special Needs:

We are proud of our proven commitment to serving the needs of our customers with the finest standing products available since 1984. We encourage you to call us with your special needs, questions, or any item that we may assist with in your standing program. Our standing and positioning specialists will be happy to answer your questions or assist you with the design and fabrication of a unit that meets your needs. Please supply as much of the information in the client measurement area of this form as possible, and list special needs in the following area to speed up the design process and / or response time.

All Prime products are made to order and are non-returnable and non-exchangeable.