



KidWalk EI

Dynamic Mobility System Patent Protected 6,832,770 B1

Canadian Order Form/Price List, effective January 2020



Client: _____ Age: _____
 Height: _____ Weight: _____ Dx: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____
 Facility: _____
 Therapist: _____ Phone: _____

PO: _____ Date: _____
 Dealer: _____
 Contact: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____ Fax: _____
 E-mail: _____

KidWalk EI

Standard Features: Auto centering pelvic guides; Lateral weight shift; Vertical weight relief with adjustable load compensating spring assembly; Angle adjustable mast; Anti-tippers; Wheel locks; Quick release wheels; Padded seat; Adjustable axle position; Light up rear casters; Two piece frame for easy transport.

KidWalk EI
\$3,905.00
Dynamic Mobility System
 Seat height 8.5"- 18"*
 Client Weight up to 65 pounds

Specifications:
 Base: 24" wide by 28" long
 Mast Height: 28"
 Upholstery: Black nylon
 Wheels: 16" Mag
 Casters: 4" Lighted



KidWalk EI additional features:
 Wrap around semi-rigid soft lateral system to fit the smaller client comfortably with a minimum of bulky hardware.
 Modular construction allows frame upgrades as the client grows.
 Ultra low seat height for true EI use and gait training.
 KidWalk EI supplies true natural gait patterns for the client.
 Minimized seat size to accommodate and support the EI client.

Frame Color:

- Candy Blue Candy Red Candy Purple Neon Pink Neon Orange Shocker Green Razzberry White Black

KidWalk Hip Support Selection: *

- Standard Configuration: Width 5" - 7" Depth 3.25" - 5.25"
 Optional Configurations: Width 5.5" - 9" Depth 3.5" - 6" 10
 Flip Up Pad Brackets: Straight (std) Degree Angle

Upper Body Support Options:

- KidWalk EI Soft Lateral System:** Flexible wraparound supports with side to side adj lateral supports with locking buckle.....\$ 682.00
Note: *Depth is variable by the wrap of the semi-rigid lateral support system, see the pictures below for reference.*
 Width: Std 4" - 7" Optional 6" - 9" with longer flex straps



EI Lateral System Open Partially closed Closed and buckled

- Rigid Lateral Chest System:** Includes lateral side supports, and swingup anterior pads adj range 6.5"-10" W x 4" - 6.5" D.....\$ 682.00
 Flip Up Chest Pad Brackets: 10 Degree Angle (std) Straight
Upper Body Soft Support Wrap: For the client that needs minimal upper body support (available in KidWalk size only).....\$ 366.30

Headrest Options:

- Headpad:** 3" x 4" flat headpad\$ 300.30
 Contoured Small Headrest: 4"H x 6"W (includes mount)\$ 458.70
 Contoured Medium Headrest: 5"H x 8"W (includes mount)\$ 531.30
 1/2 Square Tube Mount: (For use with OEM headrests)\$ 148.50
 i2i Flip Down Headrest Mount:\$ 452.10
 Specific Headrest Adapter: (Specify headrest type)..... per quote

Lower Extremity Positioning Options:

- Dynamic Leg Wraps:** Assist in supplying feedback to the client while encouraging movement in the lower extremities.....\$ 196.40
Note: *Use of the Dynamic Leg Wraps require that the seat height is set at 12 inches or above on the KidWalk EI.*
Seat Options:
Note: *All seat options are available only on a Standard KidWalk seat, the standard KidWalk EI seat cannot be modified.*
 KidWalk Seat: Replaces KidWalk EI seat\$ no charge
****Note:** *Standard seat is wider and thicker than the KidWalk EI seat, using this seat raises the seat height range to 10" - 20".***
 Wide Seat: 3" W x 8" D: Large padded seat.....\$ 270.60
 Extended Seats: Longer sides assist in controlling adduction and scissoring: Specify length of sides below\$ 422.40
 Short - 4" longer sides
 Medium - 6" longer sides
 Long - 8" longer sides

- Leg Straps:** Restricts movement of the lower extremities (Requires 6" or 8" KidWalk extended seat).....\$ 270.60
Additional Standard Seat: (When ordered with unit).....\$ 155.10
Seat size: EI Standard

Frame Options:

- Cruiser Wheels:** Allows an omni directional movement.....\$ 458.70
 Non-Lighted Caster Wheels:\$ no charge
 Sliding Ankle Prompts: Straps hold ankles to help keep the client's feet from crossing and/or to limit the client's stride\$ 357.50
 Directional Locks: Locks out caster swivel.....\$ 346.50

Wheel Options:

- Anti-Rollback Device:**.....\$ 452.10

All Prime products are made to order and are non-returnable and non-exchangeable.



Distributed in Canada by Permobil www.permobil.ca
Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944

