



KidWalk & KidWalk II

Dynamic Mobility System Patent Protected 6,832,770 B1

Canadian Order Form/Price List, effective January 2020



Client: _____ Age: _____
 Height: _____ Weight: _____ Dx: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____
 Facility: _____
 Therapist: _____ Phone: _____

PO: _____ Date: _____
 Dealer: _____
 Contact: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____ Fax: _____
 E-mail: _____

KidWalk / KidWalk II *Dynamic Mobility System*

Standard Features: Auto centering pelvic guides; Lateral weight shift; Vertical weight relief with adjustable load compensating spring assembly; Angle adjustable mast; Anti-tippers; Wheel locks; Quick release wheels; Padded seat; Adjustable axle position; Light up rear casters; Two piece frame for easy transport.

KidWalk
\$3,905.00
Dynamic Mobility System
 Seat height 12" - 22"
 Client Weight up to 65 pounds

Specifications:
 Base: 24" wide by 27" long
 Overall length: 33" - 35"
 Mast Height: 33.5"
 Upholstery: Black nylon
 Wheels: 20" Mag
 Casters: 4" Lighted



KidWalk II
\$4,724.50
Dynamic Mobility System
 Seat height 15" - 28"
 Client Weight up to 90 pounds

Specifications:
 Base: 27.5" wide by 32" long
 Overall length: 36" - 38"
 Mast Height: 37.25"
 Upholstery: Black nylon
 Wheels: 24" Mag
 Casters: 5" Lighted

Frame Color:
 Candy Blue Candy Red Candy Purple Neon Pink Neon Orange Shocker Green Razzberry White Black

KidWalk Hip Support Configuration Selection: *
 Standard Configuration: Width 5" - 9" Depth 4" x 7.5"
 Optional Configurations: Width 7" - 12" Depth 5" - 7"
 Flip Up Pad Brackets: Straight (std) 10 Degree Angle

KidWalk II Hip Support Configuration Selection: *
 Standard Configuration: Width 7"-12" Depth 5" x 7"
 Optional Width Configurations: 5.5" - 9" 8.5" - 15"
 Optional Depth Configurations: 3.5" - 6" (with non-extended hip pads)
 Flip Up Pad Brackets: Straight (std) 10 Degree Angle

Upper Body Support Options:
 Lateral Chest System: Includes lateral side supports, swingaway anterior support pads and lumbar support pad\$ **682.00**

KidWalk	KidWalk II
<input type="checkbox"/> Std 6.5" - 10.5" W x 4" - 6.5" D	<input type="checkbox"/> Std 8.5" - 15.5" W x 5.5" - 8.5" D
<input type="checkbox"/> Optional Width Adj. 8" - 14.5"	<input type="checkbox"/> Optional Width Adj. 6.5" - 10.5"
<input type="checkbox"/> Optional Depth Adj. 5.5" - 8.5"	<input type="checkbox"/> Optional Depth Adj. 4" - 6.5"

Flip Up Chest Pad Brackets: 10 Degree Angle (std) Straight

Upper Body Soft Support Wrap: For the client that needs minimal upper body support\$ **366.30**

Headrest Options:

<input type="checkbox"/> Headpad: 3" x 4" flat headpad	\$ 300.30
<input type="checkbox"/> Contoured Small Headrest: 4"H x 6"W.....	\$ 458.70
<input type="checkbox"/> Contoured Medium Headrest: 5"H x 8"W.....	\$ 531.30
<input type="checkbox"/> 1/2 Square Tube Mount: (Headrest not included).....	\$ 148.50
<input type="checkbox"/> 1/2 Flip Down Headrest Mount:	\$ 452.10
<input type="checkbox"/> i2i Flip Down Headrest Mount:	\$ 452.10

Lower Extremity Positioning Options:

<input type="checkbox"/> Extended Hip Pads: (Standard on 5" - 7" depth)	\$ 148.50
<input type="checkbox"/> Dynamic Leg Wraps: Assist in supplying feedback to the client while encouraging movement in the lower extremities.....	\$ 191.40

Seat Options:

<input type="checkbox"/> Wide Seat: 3" W x 9" D: Large padded seat.....	\$ 270.60
<input type="checkbox"/> Extended Seats: Longer sides assist in controlling adduction and scissoring: Specify length of sides below	\$ 422.40
<input type="checkbox"/> Short - 4" length	
<input type="checkbox"/> Medium - 6" length	
<input type="checkbox"/> Long - 8" length	
<input type="checkbox"/> Leg Straps: Restricts movement of the lower extremities (Requires 6" or 8" extended seat).....	\$ 270.60
<input type="checkbox"/> Additional Standard Seat: (Ordered with unit).....	\$ 155.10

Frame Options:

<input type="checkbox"/> Cruiser Wheels: Allows use as an omni directional walker for clients up to 40" tall (Not available on KidWalkII)	\$ 458.70
<input type="checkbox"/> Non-Lighted Caster Wheels:	\$ no charge
<input type="checkbox"/> Sliding Ankle Prompts: Straps hold ankles to help keep the client's feet from crossing	\$ 357.50
<input type="checkbox"/> Directional locks: Locks out caster swivel.....	\$ 346.50

Wheel Options:

<input type="checkbox"/> Anti-Rollback Device:.....	\$ 452.10
<input type="checkbox"/> 20" Spoke Wheels (KidWalk only):.....	\$ 605.00
<input type="checkbox"/> 24" Spoke Wheels (KidWalk II only):.....	\$ 605.00
<input type="checkbox"/> Spoke Protectors 20":.....	\$ 300.30
<input type="checkbox"/> Spoke Protectors 24":.....	\$ 300.30

All Prime products are made to order and are non-returnable and non-exchangeable.





Distributed in Canada by Permobil www.permobil.ca
Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944

KidWalk & KidWalk II

Dynamic Mobility System Patent Protected 6,832,770 B1

Canadian Order Form/Price List, effective January 2020

Client Name: _____



Client Measurements:

A: Chest Width ____" B: Chest Depth ____" C: Hip Width ____" D: Hip Depth ____" E: Inseam ____"

F: Top of Inseam to Axilla ____" G: Client Height ____" H: Client Weight ____lbs

Note: All Measurements Assume Symmetrical Standing Posture.

Custom Requests:

Please complete the information in the client measurement area above and list any special needs in the following area to speed up the design process and / or response time.

All Prime products are made to order and are non-returnable and non-exchangeable.