



Up-Rite I, II & III

Pediatric Mid-Line Positioning Standing Frame

Canadian Order Form/Price List, effective January 2020

Client: _____ Age: _____
 Height: _____ Weight: _____ Dx: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____
 Facility: _____
 Therapist: _____ Phone: _____

PO: _____ Date: _____
 Dealer: _____
 Contact: _____
 Address: _____
 City, Prov, PC: _____
 Phone: _____ Fax: _____
 E-mail: _____

Standard Features: Height and depth adjustable clear tabletop; Adjustable heel block; Depth, height and width adjustable knee pads; One double and one single body support strap.

UP-RITE I
\$1,424.50


Client height up to 36"

Base dimensions:
20" wide x 20" long

Tray dimensions:
18" wide x 14 1/2" deep

Knee pad cutout:
3" wide x 2 3/4" deep

Unit weight: 21lbs.



Up-Rite shown with optional caster base

UP-RITE II
\$1,853.50


Client height up to 48"

Base dimensions:
24" wide x 24" long

Tray dimensions:
21" wide x 17.5" deep

Knee pad cutout:
3" wide x 2 3/4" deep

Unit weight: 27 lbs.



Up-Rite shown with optional caster base

UP-RITE III
\$1,974.50

Client height up to 60"

Base dimensions:
32" wide x 32" long

Tray dimensions:
24" wide x 21" deep

Knee pad cutout:
3 1/2" wide x 2 3/4" deep


Unit weight: 44 lbs.



Up-Rite shown with optional caster base

- Options:**
- Wheeled caster base.....\$ **319.00**
 - Size-adjustable knee pads (pair).....\$ **284.90**
 - Double body support strap upgrade.....\$ **82.50**
 - Chest laterals (pair).....\$ **284.90**
 - Hip laterals (pair).....\$ **284.90**
 - Knee straps (pair).....\$ **126.50**

- Tray cutout with activity bowl.....\$ **214.50**
- Tray overlay Clear Opaque.....\$ **198.00**
- Foot straps (pair).....\$ **64.90**
- D-Ring foot straps (pair).....\$ **93.50**
- Adjustable head strap.....\$ **170.50**
- Shoe holders (pair, small or large).....\$ **269.50**



Client Measurements:

A: Thorax Width _____ " B: Hip Width _____ " C: Hip-to-Knee _____ "

D: Knee-to-Heel _____ " E: Shoe Size _____ " F: Shoulder Height _____ "

Note: All Measurements Assume Symmetrical Standing and Sitting Posture.

All Prime products are made to order and are non-returnable and non-exchangeable.



Distributed in Canada by Permobil www.permobil.ca
 Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944