

ROHO SOFFLEX 2 Mattress Overlay System

The SOFFLEX 2 combines cost-effectiveness with skin/soft tissue protection. The adjustable, AIR FLOATATION air cells are designed to conform to your shape, redistributing your weight across the surface. The three individual sections can be independently adjusted to address your specific support needs and your specific body and shape.

The modular mattress sections can be adjusted to a flat base or on that is anatomically shaped. The SOFFLEX 2 Mattress Overlay System is not dependent on an electrical air pump and thus it produces no sound, heat or vibrations.

The polyurethane SOFFLEX 2 Mattress Overlay is available with an enclosure cover or a reusable cover.

- ✓ Adjustable to size and weight change over time
- ✓ Assists with pain management
- ✓ Drainage holes
- ✓ Easy to clean and disinfect
- ✓ Provides 3 in. (7.5 cm) of immersion



Dimensions

Height: 3.5 in. / 9.0 cm

Width: 36 in. / 91.5 cm

Depth: 81.75 in. / 207.5 cm

Weight: 6.5 lbs. / 2.9 kg

Weight Limit

330 lbs. / 150 kg

The overlay must be properly sized to the individual and bed mattress.

Limited warranty

Mattress Overlay: 18 months

Cover: 6 months