


K300 PS Jr

permobil



Three column pricing on order/quote:

Please send order/quote to fax#: (800) 231-3256
Email to: sales@permobil.com

Dealer Information

* Contact: _____
Dealer Code: _____
* Dealer Name: _____
Address: _____
* City: _____
* Province/Postal Code: _____
* Phone#: _____
* Fax#: _____
PO#: _____
Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____
* Last Name: _____
Diagnosis: _____
Funding Source: _____
Client Age: _____

Client Measurements

* Weight: _____
* Height: _____
A) Top of Shoulders: _____
B) Chest Depth: _____
C) Chest Width: _____
D) Seat Depth: _____
E) Top of Head: _____
F) Elbow to Hand: _____
G) Seat Pan to Elbow: _____
H) Hip Width: _____
*I) Knee to Foot: _____

Permobil Ltd.
12 - 75 Mary St.
Aurora, Ontario
L4G 1G3
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.ca

Prices effective August 17, 2020 and payable in Canadian dollars.

Chair Model and Colors

Part Number	Description	MSRP - CA
<input checked="" type="checkbox"/> I103573-99-0	K300 Base PS - K300 PS0 <i>Price Includes: FWD K300 Power Base, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware for Strap Systems, Support Wheels, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 165lbs, STFH = 16.5" (17.5" with Power Tilt).</i>	8,431.00
<input type="checkbox"/> I10637	6.5mph High Speed Upgrade Package - K300 <i>Price Includes: ESP Module, High Speed 6.5mph Motor Package, and 120A Power Module.</i>	2,589.00

CHOOSE DRIVE TIRES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Tire 3x8" Grey	No Charge
<input type="checkbox"/> I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tire 3x8" Grey	No Charge

OPTIONAL DRIVE TIRES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Not available with black tires or black rims.</i>	540.00

CHOOSE SHROUD COLOR:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105701-99-0	C300/K300 Shroud - Cobalt Blue	No Charge
<input type="checkbox"/> I105702-99-0	C300/K300 Shroud - Purple Rain	No Charge
<input type="checkbox"/> I105703-99-0	C300/K300 Shroud - Onyx Black	No Charge
<input type="checkbox"/> I105704-99-0	C300/K300 Shroud - Phantom Gray	No Charge
<input type="checkbox"/> I105705-99-0	C300/K300 Shroud - Popstar Pink	No Charge
<input type="checkbox"/> I105706-99-0	C300/K300 Shroud - Volcano Red	No Charge

CHOOSE HUBCAP COLOR:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10871	Match Hubcap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i> Solid Color Hub Cap - Pair <i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Gray: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i>	No Charge

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105116-99-0	Seat Tube Fixed	No Charge
<input type="checkbox"/> I106748-99-0	Power Adjustable Seat Height-8" Travel	3,098.00

MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA
<input checked="" type="checkbox"/> I107032-99-0	Batteries, Grp 34 Gel, Pair <i>Includes Permobil Volt charger for non-Permobil gel batteries.</i>	953.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105821-99-0	Enhanced Steering Perform Unit (ESP), R-net <i>MUST be selected if ordering any ventilator tray (I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	1,223.00
<input type="checkbox"/> I105803-99-0	LED Light Kit for C/K300 w/ R-net (and New Shroud)	1,152.00

Base Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I102768-99-0	Fixed Vent Shelf for C300 <i>If selecting this part, you MUST also pick: PN I105821-99-0 - ESP Unit. Fits an LP10 Sized Vent & a U1 Battery. Must be ordered with R-net Electronics.</i>	491.00
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	540.00

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>	
<input type="checkbox"/> I108319-99-0	Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,221.00
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick - No Display <i>This Joystick does NOT have a Display Screen, does NOT contain mono jack ports, and does NOT include Bluetooth or Assignable Buttons. PC Programming required and profiles limited to five.</i>	No Charge
<input type="checkbox"/> I105700-99-0	R-net Remote Colour JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Colour LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	2,103.00

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	123.00
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	181.00
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	123.00
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	165.00
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	165.00
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	165.00
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	123.00
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	158.00
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	158.00
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	158.00

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away	404.00
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away	404.00
<input type="checkbox"/> I108037-99-0	Slimline Retractable Joystick Mount - Right	466.00
<input type="checkbox"/> I108038-99-0	Slimline Retractable Joystick Mount - Left	466.00
<input type="checkbox"/> I11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	361.00
<input type="checkbox"/> I103472-99-0	Mount Standard Joystick on Backrest Shell <i>Standard R-net Joystick is mounted to backrest. Can be chosen instead of PN ID50882. Mount on Left: _____ Mount on Right: _____</i>	323.00

Joystick Options

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Includes L-shaped joystick arm and height adjustable mounting clamp. Joystick mount sold separately (select from the mounts listed below). Please indicate if you want midline mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	711.00
<input type="checkbox"/> I108280-99-0	Bodypoint Tri-Lock Mount for R-net JS	291.00
<input type="checkbox"/> I108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick	291.00
<input type="checkbox"/> I10842	PERMOFix Mount for Mini/Compact Joystick <i>Includes a shim that allows the PERMOFix joint to be mounted on the joystick arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	291.00
Part Number	Description	MSRP - CA
<input checked="" type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with all expandable/remote joysticks.</i>	No Charge
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	No Charge
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only	1,346.00
<input type="checkbox"/> I107350-99-0	R-net Bluetooth iDevice Module <i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i>	1,087.00
<input type="checkbox"/> I107367-99-0	R-net Bluetooth Mouse Module <i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>	1,087.00
<input type="checkbox"/> I107368-99-0	Add a second Bluetooth Mouse Module <i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>	1,087.00
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,801.00

Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input checked="" type="checkbox"/> I103712-99-0	PS Junior Seat Frame <i>Weight limit 165lbs.</i>	No Charge
<input type="checkbox"/> I103715-99-0	PS Jr w/ 45° Manual Tilt <i>Must be ordered with Stroller Handles (I103811-99-0). Weight limit 165lbs.</i>	651.00
<input type="checkbox"/> I103716-99-0	PS Junior w/ 45° Power Tilt <i>Weight limit 165lbs.</i>	6,947.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I103811-99-0	Stroller Handle for PS Jr Seating System <i>Angle adjustable. This item must be selected when ordering manual tilt (I103715-99-0).</i>	225.00

Seating System Options

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below. Does not include thru drive control operation of power seat functions.</i>	No Charge
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,711.00
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>For operation of two or more seat functions through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	2,851.00

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	305.00
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	171.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,346.00
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a 315MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,284.00

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Metal Back or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10448	PS Jr Curved Metal Backrest <i>Backrest width will match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the backrest. Please note: As back height is increased, a proportional gap is created between the seat pan and the bottom of the backrest.</i>	No Charge
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge
<input type="checkbox"/> I104616-99-0	Backrest Adjustment Quick Lock, PS / PS Jr	219.00

Seat Sizing

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10432	PS Jr Seat Width 11"	No Charge
<input type="checkbox"/> I10433	PS Jr Seat Width 12"	No Charge
<input type="checkbox"/> I10434	PS Jr Seat Width 13"	No Charge
<input type="checkbox"/> I10435	PS Jr Seat Width 14"	No Charge
<input type="checkbox"/> I10436	PS Jr Seat Width 15"	No Charge
<input type="checkbox"/> I10437	PS Jr Seat Width 16"	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10985	SET SEAT DEPTH AT: _____ <i>Seat depth is adjustable from 5" to 18". Please round desired seat depth to nearest 0.5".</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10276	Miniflex Seat Cushion Leatherette <i>Approximately 2" thick. Please select cushion size: 10"W x 10"D: _____ 12"W x 12"D: _____ 14"W x 14"D: _____ Custom sizes are available upon request.</i>	607.00

Armrest Options

CHOOSE ARMREST STYLE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10824	Flip Back, Height Adjustable Armrests <i>Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (6" - 14"): _____</i>	391.00
<input type="checkbox"/> I10825	Flip Back, Height Adjustable Armrests - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (6" - 14"): _____</i>	No Charge

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I104951-99-0	3.5"x8" Unitrack Armrest Pad Right (C)	No Charge
<input type="checkbox"/> I104950-99-0	3.5"x12.5" Unitrack Armrest Pad Right (G)	No Charge
<input type="checkbox"/> I107591-99-0	4"x10" Corpus Arm w/ Permobil Gel Pad - Right	197.00
<input type="checkbox"/> I107593-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Right	203.00

CHOOSE ARMREST PAD LEFT:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I104952-99-0	3.5"x8" Unitrack Armrest Pad Left (C)	No Charge
<input type="checkbox"/> I104949-99-0	3.5"x12.5" Unitrack Armrest Pad Left (G)	No Charge
<input type="checkbox"/> I107590-99-0	4"x10" Corpus Arm w/ Permobil Gel Pad - Left	197.00
<input type="checkbox"/> I107592-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Left	203.00

Armrest Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	260.00
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	260.00
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	372.00
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	372.00
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	372.00
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	372.00
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	200.00
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	200.00
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	291.00
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	291.00

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10586	70° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate. Range of adjustment from the seat pan to the top of each footplate is: 6 - 17".</i>	No Charge
<input type="checkbox"/> I10587	90° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate. Range of adjustment from the seat pan to the top of each footplate is: 6 - 17".</i>	No Charge
<input type="checkbox"/> I10823	Aeron Legrest Assembly for PS Jr Seat - Narrow <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes 6"W x 9"D Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	404.00
<input type="checkbox"/> I11052	Aeron Legrest Assembly for PS Jr Seat - Narrow <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes 12.5"W x 10"D, Angle Adjustable, One Piece Footplate. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	404.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104684-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs</i>	No Charge
<input type="checkbox"/> I104686-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs.</i>	No Charge
<input type="checkbox"/> I107503-99-0	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	280.00

Accessories

<input type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
<input type="checkbox"/> I105849-99-0	Bodypoint Padded Hip Belt - Extra Small	148.00
<input type="checkbox"/> I104520-99-0	Bodypoint Padded Hip Belt - Medium	148.00
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	307.00
<input type="checkbox"/> I108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Med	260.00

Accessories

HEADREST OPTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I100749-99-0	MF Headrest for Koala/Ps Jr - Leatherette <i>Must also select PN I10552.</i>	508.00
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	137.00
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	258.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	137.00
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	91.00
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	91.00
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount. Note: Customization may be required if ordered along with an Omni.</i>	504.00
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. No charge with 4 power seat function configurations. Mount on left armrest: _____ Mount on right armrest: _____</i>	215.00
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	99.00
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	280.00
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	243.00
<input type="checkbox"/> I100578-99-0	Small Angle Adjustable Tray <i>Must also select PN I10556.</i>	469.00
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	137.00
<input type="checkbox"/> I108231-99-0	Permobil Adjustable Tray - 11"W	345.00
<input type="checkbox"/> I108230-99-0	Permobil Adjustable Tray - 13"W	345.00
<input type="checkbox"/> I108232-99-0	Permobil Tray Hardware Right Side Adjustable	250.00
<input type="checkbox"/> I108233-99-0	Permobil Tray Hardware Left Side Adjustable	250.00
<input type="checkbox"/> I103798-99-0	After Market Seating System Adapter <i>Used to connect after market seating systems to the UniTrack Accessory Rail on the PS Jr Seat Frame. Includes Mounting Hardware.</i>	93.00

Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106222-99-0	Stealth 3" x 5" Lateral Supports - PS Jr <i>Must also select I10723. Includes THIN Stealth hardware and specialized mounting bracket for Flip Back Armrest assembly. Only compatible with Flip Back Armrest hardware.</i>	491.00
<input type="checkbox"/> I10723	Adj Removable Stealth Lateral Hardware (Pair)	605.00

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	3,122.00
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,523.00
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,743.00

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. The Swing Away Chin Boom Assemblies are Compatible with Permobil Headrests I107470-99-0, I107469-99-0 or I107468-99-0 and must be selected with I108721-99-0.</i>	
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	296.00
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	599.00
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	959.00
<input type="checkbox"/> I10752	Haims Harness for Chin Control Mount, Large <i>For neck width up to 19cm.</i>	685.00
<input type="checkbox"/> I10753	Haims Harness for Chin Control Mount, Small <i>For neck width up to 14cm.</i>	685.00

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode & On/Off function. Switch operation force = 130 g</i> <i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i> <i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	284.00
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode & On/Off function. Switch operation force = 130 g</i>	517.00

SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	1,831.00
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	2,437.00
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i>	1,011.00
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	296.00

Total Control Head Array

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module separately.</i>	5,650.00
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module separately.</i>	No Charge
CHOOSE OCCIPITAL PAD SIZE (REQUIRED)		
<input type="checkbox"/> I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge
CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):		
<input type="checkbox"/> I10847	Fixed Mount for Occipital Pad - Head Array	No Charge
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	506.00
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):		
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	314.00
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	519.00
<input type="checkbox"/> I105925-99-0	Stealth HMO475 MultiAxis Backrest Mounting Hdwr	266.00
SELECT LEFT ARM PACKAGE (OPTIONAL)		
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	428.00
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	285.00
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	174.00
SELECT RIGHT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	428.00
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	285.00
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	174.00

Alt Drive Systems Setup

CHOOSE OMNI (REQUIRED):

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108523-99-0	Omni2 Control Module for R-net - PS Jr <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	1,995.00
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	290.00

Alt Drive Systems Setup

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	318.00
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	559.00
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	806.00
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	416.00
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	277.00
<input type="checkbox"/> I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	169.00
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	146.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I101769-99-0	Microlight Switch Qty: _____	183.00
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	138.00
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	152.00
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	135.00
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	640.00
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	1,087.00
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	39.00
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	43.00
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	47.00
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	52.00
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	43.00

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgment:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.