

TiLite Canadian Seating System Order Form

Effective August 1, 2020 (New Items in Red)

Page 1 of 5

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:		Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	Province: Code:
Ship Date:	Order #:	Phone:	Fax:

Please Email or Fax to the Permobil Manual Customer Service Team at tilitesales@permobil.com or 800.736.0925, option 2.

INSTRUCTIONS

ROHO® Seating available with new TiLite chair orders only.

SHIPPING

SSH3 Install ROHO® AGILITY® back (Must select seat back height range and set at measurement or Fixed Height when selected.).....NCO

BACK SUPPORTS

ROHO® AGILITY® ACTIVE ALUMINUM MINIMUM CONTOUR (Not available with Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

Shell Height 8" 10"

BRS14 Foam.....\$895

ROHO® AGILITY® ACTIVE ALUMINUM MID CONTOUR (Not available with Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

Shell Height 8" 10"

BRS15 Foam.....\$895

ROHO® AGILITY® MINIMUM CONTOUR (No available with Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

Shell Height 10" 13" 16" 18"

BRS5 Foam.....\$785

BRS6 Foam + Lumbar.....\$940

BRS7 Foam w/ Air Insert.....\$1,095

ACB1 AGILITY® Active Carbon Upgrade (Not available with TWIST, X Series or Bolt-on Push Handles).....\$230

BRS8 Foam w/ Air Insert + Lumbar.....\$1,095

Quick Release (Not available with AGILITY® Carbon Upgrade) OR Fixed (Not available with 18" or 20" Shell Height or Folding Chairs.).....NCO

ROHO® AGILITY® MID CONTOUR (Not available with Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

BRS9 Foam w/ Air Insert.....\$1,115

ACB2 AGILITY® Carbon Upgrade (Not available with TWIST, X Series or Bolt-on Push Handles).....\$230

BRS10 Foam w/ Air Insert + Lumbar.....\$1,115

Quick Release (Not available with AGILITY® Carbon Upgrade) OR Fixed (Not available with 18" or 20" Shell Height or Folding Chairs.).....NCO

ROHO® AGILITY® MAX CONTOUR (Not available with 10" or 13" Shell Height, Non-Folding Fixed Height Backs or Custom Height Rigidizer Bar.)

BRS11 Foam w/ Air Insert.....\$1,310

BRS12 Foam w/ Air Insert + Lumbar.....\$1,310

Quick Release OR Fixed (16" Shell Height only. Not available with Folding Chairs.).....NCO

ROHO® AGILITY™ SYSTEM ACCESSORIES

PSB1 ROHO Privacy Shield for AGILITY® Back.....\$70

COMFORT ACTA-EMBRACE™ (Not available with seat width of 13" or less.)

Shell Height 10" 12" 14" 16".....STD

SBR16 Acta-Embrace Back with Compass Hardware.....\$775

CPSB1 Comfort Privacy Shield.....\$55

COMFORT ACTA-RELIEF (Not available with seat width of 13" or less.)

Shell Height 10" 12" 14" 16".....STD

BRS13 Acta-Relief.....\$945

CPSB1 Comfort Privacy Shield.....\$55

COMFORT Install Do Not Install Install on TiLite Chair.....NCO



Permobil.com
2701 West Court Street | Pasco, Washington 99301
509.586.6117 Tel | 800.545.2266 Toll-Free
509.586.2413 Fax | 866.586.2413 Toll-Free Fax

OF0027_SeatingSystem_0419RevA

TiLite Canadian Seating System Order Form

Effective August 1, 2020 (New Items in Red)

Page 2 of 5

BACK SUPPORTS

SHORT FOLDING SEAT BACK HEIGHT RANGES AND ROHO® AGILITY™ ORIENTATION

Shell Height	Quick Release	Flipped w/ Quick Release	Fixed	Flipped w/ Fixed	Set at: (Use 1/2" increments.)
10"	<input type="checkbox"/> 16"- 18"	<input type="checkbox"/> 13"- 15"	<input type="checkbox"/> 14"- 19"	<input type="checkbox"/> 13.5"- 15.5"	Set at: _____
13"	<input type="checkbox"/> 17"- 21"	<input type="checkbox"/> 16.5"- 18.5"	<input type="checkbox"/> 17.5"- 20"	<input type="checkbox"/> 17"- 18.5"	Set at: _____
16"	<input type="checkbox"/> 19.5"- 21"	<input type="checkbox"/> 20.5"- 21.5"	<input type="checkbox"/> 19.5"- 21"	<input type="checkbox"/> 19"- 20.5"	Set at: _____
18"	<input type="checkbox"/> 22"- 24"	<input type="checkbox"/> 21.5"	N/A	N/A	Set at: _____
20"	<input type="checkbox"/> 23.5"- 25"	<input type="checkbox"/> 24.5"	N/A	N/A	Set at: _____

MEDIUM FOLDING SEAT BACK HEIGHT RANGES AND ROHO® AGILITY™ ORIENTATION

Shell Height	Quick Release	Flipped w/ Quick Release	Fixed	Flipped w/ Fixed	Set at: (Use 1/2" increments.)
10"	<input type="checkbox"/> 18"- 22"	<input type="checkbox"/> 14.5"- 18.5"	<input type="checkbox"/> 15.5"- 21"	<input type="checkbox"/> 13"- 18"	Set at: _____
AGILITY Carbon 10"	N/A	N/A	<input type="checkbox"/> 14.5"- 18"	<input type="checkbox"/> 12"- 16"	Set at: _____
13"	<input type="checkbox"/> 18"- 24.5"	<input type="checkbox"/> 16"- 21.5"	<input type="checkbox"/> 16"- 24"	<input type="checkbox"/> 16"- 21"	Set at: _____
AGILITY Carbon 13"	N/A	N/A	<input type="checkbox"/> 17"- 19"	<input type="checkbox"/> 17"- 18"	Set at: _____
16"	<input type="checkbox"/> 19"- 25.5"	<input type="checkbox"/> 19.5"- 24.5"	<input type="checkbox"/> 20"- 25"	<input type="checkbox"/> 20"- 24"	Set at: _____
AGILITY Carbon 16"	N/A	N/A	<input type="checkbox"/> 20"- 21"	<input type="checkbox"/> 19"	Set at: _____
18"	<input type="checkbox"/> 22"- 27"	<input type="checkbox"/> 22"- 25"	N/A	N/A	Set at: _____

FIXED SEAT BACK HEIGHT AND ROHO® AGILITY™ ORIENTATION (Available only with AGILITY Carbon Upgrade)

Shell Height	Fixed Back with Standard Orientation	Fixed Back with Flipped Orientation
10"	11"- 20"	13"- 20"
13"	14"- 20"	15"- 20"
16"	17"- 20"	19"- 20"

Back Orientation

Standard OR Flipped

Fixed Height

10.5" 11" 11.5" 12" 12.5" 13" 13.5" 14"
 15" 15.5" 16" 16.5" 17" 17.5" 18" 18.5"
 19" 19.5" 20"

CUSHIONS

ROHO® ADJUSTABLE SKIN PROTECTION & POSITIONING (Please see page 3 for size chart.)

CSH4 QUADTRO Select® (Fits Seat Widths 14"-19") High Profile Mid Profile Low Profile.....\$720

ROHO® ADJUSTABLE SKIN PROTECTION (Please see page 4 for size chart.)

CSH11 Single Compartment (Fits Seat Widths 14"-19") High Profile Mid Profile Low Profile.....\$720

SCD1 Sensor Ready Cushion with Smart Check Device (Not available with Low Profile.).....\$265

ROHO® AIR / FOAM COMBINATION (Please see page 3 for size chart.)

CSH14 Hybrid Elite® SR single compartment (Fits Seat Widths 14"-19").....\$745

SCD1 Smart Check Device.....\$265

CSH18 AirLITE® (Weight Limit: 275 lbs, Fits chair widths 13"-20").....\$250

COMFORT CUSHIONS

CSH2 Comfort Curve Cushion (Stretch-Air Cover. Not available with Seat Taper).....\$155

CSH21 Versa™ X (Not available with Seat Taper. Please see below for size chart.).....\$615

CSH23 Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion (Please see below for size chart.).....\$470

CSH24 Comfort Elements™ 2" Thick Comfort-Tek Cushion (Limited width and depth sizes. Please see below for size chart.).....\$140

TiLite Seating System Order Form

Effective August 1, 2020 (New Items in Red)

Page 3 of 4

CUSHIONS

Versa™ X and Acta-Embrace

Rear Seat Width	Seat Depth
12"	10", 11", 12", 13", 14", 15", 16"
13"	10", 11", 12", 13", 14", 15", 16", 17"
14"	10", 11", 12", 13", 14", 15", 16", 17", 18"
15"	14", 15", 16", 17", 18", 19", 20"
16"	14", 15", 16", 17", 18", 19", 20"
17"	15", 16", 17", 18", 19", 20"
18"	15", 16", 17", 18", 19", 20"
19"	16", 17", 18", 19", 20"
20"	16", 17", 18", 19", 20"

Hybrid Elite®

Rear Seat Width	Seat Depth
14"	14", 15", 16"
15"	15", 16", 17"
16"	15", 16", 17", 18", 20"
17"	16", 17", 18"
18"	16", 17", 18", 20"

AirLITE®

Rear Seat Width	Seat Depth
13"	13", 15"
14"	14", 16"
15"	15", 17"
16"	16", 18"
17"	17", 18", 20"
18"	16", 18", 20"
19"	17", 19"
20"	17", 18"

Elements Gel Cushion

Rear Seat Width	Seat Depth
16"	16", 18", 20"
18"	16", 18", 20"
20"	16", 18", 20"

QUADTRO SELECT®

Rear Seat Width	12" Seat Depth	13" Seat Depth	14" Seat Depth	15" Seat Depth	16" Seat Depth	17" Seat Depth	18" Seat Depth	19" Seat Depth	20" Seat Depth
12"	Low, High	Low, High	N/A	N/A	N/A	N/A	N/A	N/A	N/A
13"	N/A	Low, High	Low, High	Low, High	Low, High	Low, High	N/A	N/A	N/A
14"	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High	Low, High	N/A
15"	N/A	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High	Low, High	N/A
16"	N/A	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High
17"	N/A	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High
18"	N/A	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High
19"	N/A	N/A	Low, High	Low, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, Mid, High	Low, High

TiLite Seating System Order Form

Effective August 1, 2020 (New Items in Red)

CUSHIONS

Single Compartment

Italicized profiles are available without ROHO Smart Check

Bolded profiles are available with ROHO Smart Check

Rear Seat Width	12" Seat Depth	13" Seat Depth	14" Seat Depth	15" Seat Depth	16" Seat Depth	17" Seat Depth	18" Seat Depth	19" Seat Depth	20" Seat Depth
12"	Low, High	Low, High	N/A	N/A	N/A	N/A	N/A	N/A	N/A
13"	N/A	<i>Low, Mid, High</i> High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, High</i> Mid, High	<i>Low, High</i> Mid, High	<i>Low, High</i> High	N/A	N/A
14"	N/A	<i>Low, High</i> High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, High</i> High	<i>Low, High</i>	N/A
15"	N/A	<i>Low, High</i> High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, High</i> High	<i>Low, High</i>	N/A
16"	N/A	N/A	<i>Low, High</i> High	<i>Low, High</i> High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, High</i>
17"	N/A	N/A	<i>Low, High</i> High	<i>Low, High</i> High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, High</i>
18"	N/A	N/A	<i>Low, High</i> Mid, High	<i>Low, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High
19"	N/A	N/A	<i>Low, High</i> Mid, High	<i>Low, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High
20"	N/A	N/A	N/A	N/A	<i>Low, High</i> High	<i>Low, High</i> High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High	<i>Low, Mid, High</i> Mid, High

GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Bodypoint® is a registered trademark of Bodypoint Inc. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc.. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Strain® is a registered trademark of Q'Strain Inc.

NOTES
