

Hälsa bed



Ideal Home Care Bed

The Hälsa bed from Permobil is an electric, adjustable home care bed that comes in regular and plus sizes. The Trendelenburg and reverse Trendelenburg feature comes standard. Additionally, it offers the capability of moving both the head and foot sections simultaneously, an upgraded caster option, and easy disassembly for transportation and storage. Hälsa regular bed has removable composite deck sections for easy cleaning.

Standard features include:

- Dual posts on each side increases stability at all heights
- Low Position: 9.875" (12.125" with large caster)
- High Position: 25.625" (27.875" with large caster)
- Trendelenburg/reverse Trendelenburg
- Virtually noise free with low electrical consumption DC Linak motors
- Rolls at any height
- Hooks for incontinence bag
- 3 locking casters & 1 directional locking caster
- Side and foot mattress retainers
- Receptacle for trapeze bar
- 2 receptacles for IV holder stems
- Disassembles into 4 sections for transport/storage

At a Glance

Regular size (36" x 80") and Plus size (4 width expansion options: 36", 39", 42" & 48")

Low height of 9.875", high height of 25.625"

Weight capacity: Hälsa - 375 lbs, Hälsa Plus - 500 lbs

Trendelenburg and reverse Trendelenburg feature comes standard



Optional Features

- Half-side rails - interchangeable between head and foot sections
- Full side rails
- Assist Bar - interchangeable between head and foot sections
- Trapeze bar
- Nature and Wine colour Head/Footboards
- Battery backup
- Upgraded large caster option

Specifications

Size	Available in 36" x 80", and Plus size with 4 width expansion options: 36", 39", 42" & 48"
Product Weight	265 lbs
Patient Weight Limit	Reg. is 375 lbs; Plus is 500 lbs
Warranties	Frame (10 years), Mechanical (4 years), Electrical (2 years), Wood Products (1 year)