

Tröst bed



At a Glance

Low height of 6.75" and a high height of 30"

Weight capacity of 600 lbs

Centre step safety lock at foot end

12-button hand pendant with Trendelenburg, reverse Trendelenburg and Comfort Chair



Ideal Heights for Care

The Tröst bed has a low height of 6.75" and a high height of 30" to meet all your caregiving and client needs. With a 12-button hand pendant and optional footboard attendant controls, customized comfort is an easy adjustment.

Standard features include:

- 12-button hand pendant with Trendelenburg, reverse Trendelenburg and Comfort Chair
- Centre step safety lock at foot end
- Flip down directional caster locks at head end
- Quick release Head/Footboard hardware
- Swivel locking casters
- Ratchet foot section
- Wall bumpers
- Mattress retainers

Fixed Freight Policy

With one-person delivery, Tröst beds are simple to ship and set up. Whether you buy one bed or over 10, Permobil has a special fixed freight policy - simply call for a quote.

Want to save even more money? We offer free freight on any mattress or support surface ordered with a bed.

Optional Features

- Tröst Plus model which is expandable to 39" and 42"
- Length extension from 80" to 84"
- Embedded footboard attendant controls
- Half-side rails - interchangeable between head and foot sections
- Assist Bar - interchangeable between head and foot sections
- Trapeze
- Nature and Wine colour Head/Footboards
- Transport trolley
- Battery backup



Specifications

Size	Available in 35", 39" and 42" width 80" or 84" lengths offered
Product Weight	250 lbs
Patient Weight Limit	600 lbs
Warranties	Frame (10 years), Mechanical (4 years), Electrical (2 years), Wood Products (1 year)