

# M5 Corpus

# permobil



## Dealer Information

\* **Contact:** \_\_\_\_\_  
Dealer Code: \_\_\_\_\_  
\* **Dealer Name:** \_\_\_\_\_  
Address: \_\_\_\_\_  
\* **City:** \_\_\_\_\_  
\* **Province/Postal Code:** \_\_\_\_\_  
\* **Phone#:** \_\_\_\_\_  
\* **Fax#:** \_\_\_\_\_  
PO#: \_\_\_\_\_  
Email Address: \_\_\_\_\_


## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* **First Name:** \_\_\_\_\_  
\* **Last Name:** \_\_\_\_\_  
Diagnosis: \_\_\_\_\_  
Funding Source: \_\_\_\_\_  
Client Age: \_\_\_\_\_

## Client Measurements

\* **Weight:** \_\_\_\_\_  
\* **Height:** \_\_\_\_\_  
A) Top of Shoulders: \_\_\_\_\_  
B) Chest Depth: \_\_\_\_\_  
C) Chest Width: \_\_\_\_\_  
D) Seat Depth: \_\_\_\_\_  
E) Top of Head: \_\_\_\_\_  
F) Elbow to Hand: \_\_\_\_\_  
G) Seat Pan to Elbow: \_\_\_\_\_  
H) Hip Width: \_\_\_\_\_  
\* **I) Knee to Foot:** \_\_\_\_\_



Three column pricing on order/quote:

Please send order/quote to fax#: (800) 231-3256  
Email to: [sales@permobil.com](mailto:sales@permobil.com)

**Permobil Ltd.**  
12 - 75 Mary St.  
Aurora, Ontario  
L4G 1G3  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.ca](http://www.permobil.ca)

Prices effective **January 1, 2022** and payable in Canadian dollars.

## Chair Model and Colors

Part Number	Description	MSRP - CA
<input type="checkbox"/> I11353	M5 Base Corpus - M5 SPO <i>Price Includes: MWD M5 Power Base Set Up For One Power Seat Function, Corpus Seat Frame, FlexLink Pro Suspension, 8" Flat-Free Casters, Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 7.5mph, Weight Capacity = 300lbs.</i>	9,821.00
<input type="checkbox"/> I109053-99-0	M5 Base Corpus - M5 MPO <i>Price Includes: MWD M5 Power Base Set Up For Multiple Power Seat Functions, Corpus Seat Frame, FlexLink Pro Suspension, 8" Flat-Free Casters, Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 7.5mph, Weight Capacity = 300lbs.</i>	9,821.00
<input type="checkbox"/> I108852-99-0	Permobil Connect <i>By selecting, user agrees to activate the Permobil Connect-ready wheelchair for connected services.</i>	No Charge

## CHOOSE DRIVE TIRES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108414-99-0	Flat Free Tire 3"W x 14"T Black - Standard Tread <i>Designed for a smooth ride on indoor and outdoor surfaces. Includes flat free standard tread caster wheels.</i>	No Charge
<input type="checkbox"/> I108411-99-0	Pneumatic Tire 3"W x 14"T Black - Standard Tread <i>Includes flat free standard tread caster wheels.</i>	No Charge
<input type="checkbox"/> I108413-99-0	Flat Free Tire 3"W x 14"T Black - Aggressive Tread <i>Includes flat free aggressive tread caster wheels.</i>	No Charge

## OPTIONAL DRIVE TIRES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108412-99-0	Pneumatic/Snow Tire 3"W x 14"T Black - Aggr Tread <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>	567.00

## CHOOSE SHROUD/RIM COLOR:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I109117-99-0	M5 Color Kit - Aqua	No Charge
<input type="checkbox"/> I109118-99-0	M5 Color Kit - Sky	No Charge
<input type="checkbox"/> I109119-99-0	M5 Color Kit - Violet	No Charge
<input type="checkbox"/> I109120-99-0	M5 Color Kit - Flamingo	No Charge
<input type="checkbox"/> I109121-99-0	M5 Color Kit - Chili	No Charge
<input type="checkbox"/> I109122-99-0	M5 Color Kit - Canary	No Charge
<input type="checkbox"/> I109123-99-0	M5 Color Kit - Spring	No Charge
<input type="checkbox"/> I109124-99-0	M5 Color Kit - Arctic	No Charge
<input type="checkbox"/> I109125-99-0	M5 Color Kit - Sunset	No Charge
<input type="checkbox"/> I109177-99-0	M5 Color Kit - Midnight	No Charge
<input type="checkbox"/> I109126-99-0	M5 Color Kit - Gold	No Charge
<input type="checkbox"/> I109127-99-0	M5 Color Kit - Bronze	No Charge
<input type="checkbox"/> I109128-99-0	M5 Color Kit - Platinum	No Charge
<input type="checkbox"/> I109130-99-0	M5 Color Kit - Carbon Fiber	No Charge
<input type="checkbox"/> I109131-99-0	M5 Color Kit - Mossy Oak Break-Up Country	No Charge
<input type="checkbox"/> I109132-99-0	M5 Color Kit - CDG Poly	No Charge
<input type="checkbox"/> I109134-99-0	M5 Color Kit - Fireball	No Charge
<input type="checkbox"/> I109135-99-0	M5 Color Kit - Water Lily	No Charge

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
-------------	-------------	-----------

### Chair Model and Colors

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	Hub Cap Color: <i>Please write in Hub Cap Color below. Hub Cap Color will be matched to the Shroud Color if no selection is made.</i> _____	

### Base Options

#### MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107086-99-0	Active Height Power Adj. Seat Height - 14" Travel <i>Must also select Power Tilt when selecting this option.</i>	3,721.00
<input type="checkbox"/> I11117	Seat Tube Fixed	No Charge

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107435-99-0	Set Seat Height to 17.5"	No Charge
<input type="checkbox"/> I107433-99-0	Set Seat Height at 18.5" <i>If ordering Power Adjustable Seat Height, 1" of elevator travel is used to set seat at 18.5".</i>	No Charge
<input type="checkbox"/> I107434-99-0	Set Seat Height at 19.5" <i>If ordering Power Adjustable Seat Height, 2" of elevator travel is used to set seat at 19.5".</i>	No Charge

#### MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA
<input checked="" type="checkbox"/> I108999-99-0	BATTERY OPTIONS: Permobil Batteries, Grp 24 AGM GEL, Pair <i>Includes Permobil VoltPro charger optimized for Permobil batteries.</i>	1,013.00

#### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107998-99-0	Enhanced Steering Perform Unit (ESP), R-net <i>Provided at no charge when ordered with an Omni2 Control Module.</i>	1,284.00
<input type="checkbox"/> I109170-99-0	Essentials Box with Reflectors <i>Please Note: The total backrest angle will be limited to 170°, or further depending on other backrest-mounted items, to prevent contact with the Essentials Box.</i>	580.00
<input type="checkbox"/> I109169-99-0	Essentials Box with Lights <i>Please Note: The total backrest angle will be limited to 170°, or further depending on other backrest-mounted items, to prevent contact with the Essentials Box. Must also select Light Kit.</i>	965.00
<input type="checkbox"/> I109168-99-0	LED Light Kit for M3 / M5 <i>No charge with 4 power seat function configurations.</i>	1,255.00
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 330mm Tall, adjustable Width between 200mm - 330mm, and adjustable Depth between 50mm - 165mm. Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option. Additional drive restrictions may apply.</i>	1,819.00
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray <i>Mounts to rear of Laptop Ventilator Tray. Measures 290mm W x 170mm H x 42mm D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	291.00
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 175mm Tall. Width infinitely adjustable 90mm - 170mm. Depth infinitely adjustable 95mm - 125mm. Custom Option.</i>	411.00
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	567.00

### Joystick Options

#### CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA
-------------	-------------	-----------

## Joystick Options

## CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>	
<input type="checkbox"/> I107779-99-0	Color Joystick Module, (CJSM) <i>This R-net color joystick features: A Large High-Resolution 3.5" VGA color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,444.00
<input type="checkbox"/> I107778-99-0	Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,282.00
<input type="checkbox"/> I11186	Permobil Joystick Module, NC <i>Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display.</i>	No Charge
<input type="checkbox"/> I107829-99-0	3.5mm Y-Splitter Cable for PJSM/CJSM <i>For use with the PJSM or CJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	45.00
<input type="checkbox"/> I105700-99-0	R-net Remote Colour JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Colour LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	2,208.00

## CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	129.00
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	190.00
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	129.00
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	173.00
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	173.00
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	173.00
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	129.00
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	166.00
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	166.00
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	166.00

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10664	Joystick Mount Right - Fixed	No Charge
<input type="checkbox"/> I10665	Joystick Mount Left - Fixed	No Charge
<input type="checkbox"/> I108037-99-0	Slimline Retractable Joystick Mount - Right	489.00
<input type="checkbox"/> I108038-99-0	Slimline Retractable Joystick Mount - Left	489.00
<input type="checkbox"/> I11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	379.00
<input type="checkbox"/> I105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN I105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	339.00
<input type="checkbox"/> I108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Includes L-shaped joystick arm and height adjustable mounting clamp. Joystick mount sold separately (select from the mounts listed below). Please indicate if you want midline mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	747.00
<input type="checkbox"/> I108280-99-0	Bodypoint Tri-Lock Mount for R-net JS	306.00
<input type="checkbox"/> I108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick	306.00

## Joystick Options

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10842	PERMOFix Mount for Mini/Compact Joystick <i>Includes a shim that allows the PERMOFix joint to be mounted on the joystick arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	306.00
<input checked="" type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with all expandable/remote joysticks.</i>	No Charge
<input checked="" type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	No Charge
<input type="checkbox"/> I105225-99-0	Co-Pilot Attendant Control R-net - Corpus	5,444.00
<input type="checkbox"/> I105408-99-0	Attendant Control R-net	1,413.00
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,891.00

## Seating System

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>Speed limitations are imposed when the total backrest angle is greater than 150° or when the seat is tilted forward more than 5°. Drive is inhibited at 40° of tilt for 17.5" seat heights and user weights &lt;220lbs. For taller seat heights and user weights &gt;220lbs drive inhibit occurs at 35° of tilt.</i>	
<input type="checkbox"/> I11211	Corpus Seat w/ Manual Recline <i>Manual recline angle limited to 120°.</i>	2,620.00
<input type="checkbox"/> I11212	Corpus Seat w/ No Tilt & 180° Power Recline	9,908.00
<input type="checkbox"/> I11213	Corpus Seat w/ 50° Tilt	8,884.00
<input type="checkbox"/> I11214	Corpus Seat w/ 50° Power Tilt & 180° Power Recline	14,638.00

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>Active Reach should be recommended by a licensed clinical professional and/or ATP. While operating Active Reach, positioning equipment such as lap belts and chest supports should be considered to meet the user's specific safety needs. Note: A low speed limitation is imposed from 5° to 10° of Active Reach, an extra low speed limitation is imposed from 10° to 20° of Active Reach, and a drive lockout is imposed above 20° of Active Reach.</i>	
<input type="checkbox"/> I11141	5° Active Reach (Anterior Tilt) <i>When ordered with a Power Adjustable Seat Height, tilt range will be -5° to +50°. When ordered without a Power Adjustable Seat Height, tilt range will be -5° to +45°. If this option is selected with a fixed seat, the minimum seat to floor height at 0° of tilt will be approximately 19.5". Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 97° when manual legrest is selected.</i>	No Charge
<input type="checkbox"/> I107090-99-0	10° Active Reach (Anterior Tilt) <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height and Power Tilt. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 103° when manual legrest is selected.</i>	No Charge
<input type="checkbox"/> I107091-99-0	20° Active Reach (Anterior Tilt) <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	3,199.00



## Seating System Options

### CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Does not include thru drive control operation of power seat functions.</i>	No Charge
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	1,797.00
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>For operation of two or more seat functions through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	2,994.00

### CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA
PLEASE NOTE: <i>ICS switchbox is only available with chairs that include standing, 20° or more of ActiveReach or an R-Net LED Remote Joystick.</i>		
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	320.00
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	180.00

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,413.00
<input type="checkbox"/> I108416-99-0	Remote Stop System <i>Consists of a 315MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>	1,556.00
<input type="checkbox"/> I10750	Memory Seat Program, Corpus w/ R-net	131.00
<input type="checkbox"/> I106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	53.00
<input type="checkbox"/> I10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge

## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge
<input type="checkbox"/> I105066-99-0	Fixed Backrest Kit - No Mechanical Shear	58.00

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
-------------	-------------	-----------

## Backrest Options

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	<b>BACKREST HEIGHT NOTE:</b> <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Corpus Ergo Backs include lateral wedges and foam lumbar pads. Solid Backs do not include foam lumbar pads.</i>	
<input type="checkbox"/> I106545-99-0	Corpus Ergo Back - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,159.00
<input type="checkbox"/> I106546-99-0	Corpus Ergo Back - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,159.00
<input type="checkbox"/> I106547-99-0	Corpus Ergo Back - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,159.00
<input type="checkbox"/> I106548-99-0	Corpus Ergo Back - 20"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,159.00

## PLEASE NOTE:

*Taller backrests listed below include: lateral wedges, foam lumbar pads\*, UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided. \* Foam lumbar pads not included with Solid Backs.*

<input type="checkbox"/> I106532-99-0	Corpus Ergo Back - 14"W x 23" - 28"T (C3GERGBCA)	1,159.00
<input type="checkbox"/> I106533-99-0	Corpus Ergo Back - 16"W x 23" - 28"T (C3GERGBCA)	1,159.00
<input type="checkbox"/> I106534-99-0	Corpus Ergo Back - 18"W x 23" - 28"T (C3GERGBCA)	1,159.00
<input type="checkbox"/> I106535-99-0	Corpus Ergo Back - 20"W x 23" - 28"T (C3GERGBCA)	1,159.00

## MUST ALSO SELECT BACK HEIGHT:

Set Corpus Backrest at 23" \_\_\_\_\_  
 Set Corpus Backrest at 24" \_\_\_\_\_  
 Set Corpus Backrest at 25" \_\_\_\_\_  
 Set Corpus Backrest at 26" \_\_\_\_\_  
 Set Corpus Backrest at 27" \_\_\_\_\_  
 Set Corpus Backrest at 28" \_\_\_\_\_

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I11335	Leatherette Cushion - Corpus Ergo Back <i>Includes Dual Density Backrest Foam.</i>	No Charge
<input type="checkbox"/> I11336	Stretch-Air Cushion - Corpus Ergo Back <i>Includes Dual Density Backrest Foam with Stretch-Air breathable fabric cover.</i>	No Charge

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	<b>ROHO AGILITY NOTE:</b> <i>The ROHO AGILITY Back Systems are available in widths of 14", 16", 18" and 20" and shell heights of 16", 18" and 20". Direct Mount system includes backrest shell, cushion, air inserts, UT rail for headrest and Direct Mount adapter for Corpus.</i>	
<input type="checkbox"/> I11258	ROHO AGILITY Max Contour Back System - 16" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 22" and 25". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,167.00
<input type="checkbox"/> I11259	ROHO AGILITY Max Contour Direct Mount - 18" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 22" and 26". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,167.00
<input type="checkbox"/> I11260	ROHO AGILITY Max Contour Direct Mount - 20" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 24" and 27". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,167.00

## MUST ALSO SELECT BACK HEIGHT:

Set AGILITY Backrest at 22" \_\_\_\_\_  
 Set AGILITY Backrest at 23" \_\_\_\_\_

## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	Set AGILITY Backrest at 24" _____	
	Set AGILITY Backrest at 25" _____	
	Set AGILITY Backrest at 26" _____	
	Set AGILITY Backrest at 27" _____	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I109029-99-0	BodiLink Premium 3" x 5" Lateral Supports - Pair <i>Must also select I109031-99-0. Includes foam wedge for additional support. Only compatible with Corpus Ergo backrest shells.</i>	470.00
<input type="checkbox"/> I109030-99-0	BodiLink Premium 4" x 6" Lateral Supports - Pair <i>Must also select I109031-99-0. Includes foam wedge for additional support. Only compatible with Corpus Ergo backrest shells.</i>	470.00
<input type="checkbox"/> I109031-99-0	Swing-Away Adj BodiLink Lateral Hdwr Corpus Pair	725.00

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	516.00
<input type="checkbox"/> I106072-99-0	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrix Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	516.00
<input type="checkbox"/> I10725	Split Backrest Cushion - Corpus <i>Custom Option.</i>	384.00
<input type="checkbox"/> I10786	Tailored Backrest Cushion - Corpus <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion when ordering this custom option.</i>	120.00

## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107023-99-0	Set Seat Width at 17": _____	No Charge
<input type="checkbox"/> I107024-99-0	Set Seat Width at 19": _____	No Charge
<input type="checkbox"/> I107025-99-0	Set Seat Width at 21": _____	No Charge
<input type="checkbox"/> I107026-99-0	Set Seat Width at 23": _____	No Charge

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>There are three seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>	
<input type="checkbox"/> I107013-99-0	Seat Depth 14"	No Charge
<input type="checkbox"/> I107014-99-0	Seat Depth 15"	No Charge
<input type="checkbox"/> I107015-99-0	Seat Depth 16"	No Charge
<input type="checkbox"/> I107016-99-0	Seat Depth 17"	No Charge
<input type="checkbox"/> I107017-99-0	Seat Depth 18"	No Charge
<input type="checkbox"/> I107018-99-0	Seat Depth 19"	No Charge
<input type="checkbox"/> I107019-99-0	Seat Depth 20"	No Charge
<input type="checkbox"/> I107020-99-0	Seat Depth 21"	No Charge



## Seat Sizing

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107021-99-0	Seat Depth 22"	No Charge

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
NOTE: A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.		
<input type="checkbox"/> I11337	Corpus Ergo Cushion - Leatherette Cover	654.00
<input type="checkbox"/> I11338	Corpus Ergo Cushion - Stretch-Air Cover	654.00
<input type="checkbox"/> I10965	Corpus Ergo Air Cushion (C3GRHOSCA) <i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23". Not Available in 21" x 17"</i>	694.00
<input type="checkbox"/> I11113	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part. Part Number will be provided on the quote / order confirmation based on the Seat Width and Seat Depth selected.</i>	208.00
<input type="checkbox"/> I11166	ROHO QUADTRO SELECT HIGH PROFILE	777.00
<input type="checkbox"/> I11168	ROHO HIGH PROFILE Single Compartment	636.00
<input type="checkbox"/> I11170	ROHO HIGH PROFILE Sensor Ready <i>Single compartment cushion designed to work with Smart Check Device I107505-99-0. Not available in all seat widths and depths.</i>	636.00
<input type="checkbox"/> I11173	Incontinence Cover for ROHO Cushion <i>Cover size will be matched to ROHO cushion size. Sent as a part.</i>	145.00
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge

## Armrest Options

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
PLEASE NOTE: <i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>		
<input type="checkbox"/> I106543-99-0	Narrow Armbar Kit, Inf Adjustable 15"W - 19"W	No Charge
<input type="checkbox"/> I106544-99-0	Wide Armrest Kit, Inf Adjustable 19"W - 23"W <i>Required when selecting a 20" wide backrest.</i> Set the Width Between the Armrests to: _____ <i>Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>	No Charge

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105069-99-0	Height Adjustable Armrests - NC SET ARMREST HEIGHT AT: _____ <i>Corpus Armrests are infinitely adjustable from 8" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>	No Charge

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105071-99-0	Left Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge
<input type="checkbox"/> I105070-99-0	Right Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge

## Armrest Options

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge
<input type="checkbox"/> I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105105-99-0	4"x16" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105999-99-0	4"x18" Corpus HD Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I107214-99-0	4"x16" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I105995-99-0	4"x18" Corpus HD Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I107593-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Right	213.00
<input type="checkbox"/> I107595-99-0	4"x16" Corpus Arm w/ Permobil Gel Pad - Right	218.00
<input type="checkbox"/> I107597-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Right	227.00
<input type="checkbox"/> I108861-99-0	Moldable Forearm Support - Right 11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.	360.00
<input type="checkbox"/> I108863-99-0	Contoured Hand Pad Medium - Right 6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.	179.00
<input type="checkbox"/> I108865-99-0	Contoured Hand Pad Large - Right 7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.	189.00

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105107-99-0	4"x16" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105998-99-0	4"x18" Corpus HD Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I107215-99-0	4"x16" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I105994-99-0	4"x18" Corpus HD Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I107592-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Left	213.00
<input type="checkbox"/> I107594-99-0	4"x16" Corpus Arm w/ Permobil Gel Pad - Left	218.00
<input type="checkbox"/> I107596-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Left	227.00
<input type="checkbox"/> I108860-99-0	Moldable Forearm Support - Left 11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.	360.00
<input type="checkbox"/> I108862-99-0	Contoured Hand Pad Medium - Left 6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.	179.00
<input type="checkbox"/> I108864-99-0	Contoured Hand Pad Large - Left 7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.	189.00

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105463-99-0	Corpus Arm Trough Adapter 1" - Pair For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.	375.00
<input type="checkbox"/> I105462-99-0	Corpus Arm Trough Adpater 7/8" - Pair For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option. Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	375.00

### Armrest Options

#### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	273.00
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	273.00
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	391.00
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	391.00
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	391.00
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	391.00
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	210.00
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	210.00
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	306.00
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	306.00

### Legrest Options

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
CHOOSE LEGREST TYPE:		
<input type="checkbox"/> I107763-99-0	Corpus Manual Elevating Legrest	No Charge
<input type="checkbox"/> I107926-99-0	Power Elevating Legrest - 170°	3,838.00
CHOOSE FOOTPLATE STYLE:		
<input type="checkbox"/> I11197	2pc Narrow Corpus Footplates 9"D x 6"W	No Charge
<input type="checkbox"/> I11196	2pc Wide Corpus Footplates - 9"D x 7.5"W	No Charge
<input type="checkbox"/> I106760-99-0	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	410.00
CALF SUPPORT OPTIONS:		
<input type="checkbox"/> I105252-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge
<input type="checkbox"/> I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	424.00
<input type="checkbox"/> I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad.</i>	424.00

#### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
CORPUS VS LEGREST AND KNEE SUPPORT:		
<i>NOTE: The Corpus VS Legrest is only available with 20°, 30° and 45° of Active Reach.</i>		
<input type="checkbox"/> I107925-99-0	Corpus VS Power Elevating Legrest - 180° <i>Includes power articulation.</i>	3,942.00
<input type="checkbox"/> I107038-99-0	Knee Support - Corpus VS <i>Select when ordering Corpus VS Legrest. Must also select I11129.</i>	226.00

## Legrest Options

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I11129	Adj Removable Knee Support Hardware (1 Each)	330.00
CHOOSE FOOTPLATE STYLE:		
<input type="checkbox"/> I11198	2pc Narrow Corpus VS Footplates - 9"D x 6"W	No Charge
CALF SUPPORT OPTIONS:		
<input type="checkbox"/> I107211-99-0	VS Calf Support Kit (J) w Mesh Cover & Adj Hardware 5"W x 6.5"T Pad. Compatible with footplate heights of 10" - 18".	No Charge
<input type="checkbox"/> I106190-99-0	VS Calf Support Kit (A) w Mesh Cover & Adj Hardware 5.5"W x 6"T Pad. Compatible with footplate heights of 10" - 18".	No Charge
<input type="checkbox"/> I106188-99-0	VS Calf Support Kit (H) w Mesh Cover & Adj Hardware 7"W x 7"T Curved Pad. Compatible with footplate heights of 11" - 18".	455.00
<input type="checkbox"/> I106186-99-0	VS Calf Support Kit (D) w Mesh Cover & Adj Hardware 8"W x 8"T Pad. Compatible with footplate heights of 12" - 18".	455.00

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9".</i>	

## CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106530-99-0	Stealth Gel Pads - Narrow Footplates (Pair)	432.00
<input type="checkbox"/> I106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	444.00

## Accessories

<input type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
<input type="checkbox"/> I104520-99-0	Bodypoint Padded Hip Belt - Medium	155.00
<input type="checkbox"/> I104514-99-0	Bodypoint Non-Padded Hip Belt - Large	125.00
<input type="checkbox"/> I104516-99-0	Bodypoint Padded Hip Belt - Large	213.00
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	322.00
<input type="checkbox"/> I100600-99-0	Positioning Belt Retractable Long	350.00
<input type="checkbox"/> I108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Med	273.00
<input type="checkbox"/> I108266-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Large	289.00

## HEADREST OPTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107470-99-0	Permobil 7"W Headrest <i>Must also select headrest mounting hardware.</i>	295.00
<input type="checkbox"/> I107469-99-0	Permobil 10"W Headrest <i>Must also select headrest mounting hardware.</i>	295.00
<input type="checkbox"/> I107468-99-0	Permobil 14"W Headrest <i>Must also select headrest mounting hardware.</i>	295.00
<input type="checkbox"/> I107480-99-0	Adj/Removable Headrest Link Hardware - Standard	315.00
<input type="checkbox"/> I109178-99-0	Extra Link for Link Headrest Hardware <i>Used to extend range of standard link hardware. Sent as part.</i>	63.00
<input type="checkbox"/> I108721-99-0	Adj/Removable BodiLink Head Support Hdwr - MAX 8 <i>Must be selected when ordering more than 2 swing-away headrest accessories or chin boom assemblies.</i>	325.00
<input type="checkbox"/> I108722-99-0	BodiLink Headrest Bracket - Corpus 23-28 <i>For 23" - 28" Corpus Ergo backrest heights.</i>	148.00

## Accessories

## HEADREST OPTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108723-99-0	Bodilink Headrest Bracket Corpus Short/AGILITY <i>For 20" Corpus Ergo backrests and AGILITY Max Contour Back Systems.</i>	148.00
<input type="checkbox"/> I11305	Swing Away Egg Switch Assembly <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	645.00
<input type="checkbox"/> I11306	Swing Away Microlight Assembly <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	740.00
<input type="checkbox"/> I11307	Swing Away Facial Pad Assembly 4"L x 2"T <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	525.00
<input type="checkbox"/> I11308	Swing Away Facial Pad Assembly 6"L x 2"T <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	525.00
<input type="checkbox"/> I11276	Provide Headrest Receiver for AGILITY Backrest <i>Select for AGILITY backs and third party backs ordered off chair.</i>	No Charge
<input type="checkbox"/> I106331-99-0	Permobil 10"W Headrest, Leatherette <i>Not available with ROHO AGILITY Back Systems. Must also select PN I10552.</i>	672.00
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	144.00
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	271.00

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	403.00
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	403.00
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left:_____ Mount Amp Suppt on Right:_____</i>	144.00
<input type="checkbox"/> I109040-99-0	BodiLink Premium LPTS Pad (B), Stretch-Air - Pair <i>4.5"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	284.00
<input type="checkbox"/> I109043-99-0	BodiLink Premium LPTS Pad (A), Stretch-Air - Pair <i>6"W x 5.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	284.00
<input type="checkbox"/> I109041-99-0	BodiLink Premium LPTS Pad (C), Stretch-Air - Pair <i>8"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	284.00
<input type="checkbox"/> I109042-99-0	BodiLink Premium LPTS Pad (G), Stretch-Air - Pair <i>12"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	284.00
<input type="checkbox"/> I109044-99-0	Adj/Removable BodiLink LPTS Hdwr Med Corpus Pair <i>Medium TT Power Mount hardware can be extended upward to clear cushions of 3-5.5" thick.</i>	336.00
<input type="checkbox"/> I109194-99-0	Adj/Removable Lateral Thigh/Hip Hardware - Pair	144.00
<input type="checkbox"/> I106762-99-0	Abductor Kit Corpus - Leatherette <i>Not compatible with Corpus VS Legrest. Must also select PN I10555.</i>	399.00
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	144.00
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus Seat	554.00
<input type="checkbox"/> I105292-99-0	Push Handles, Corpus Seat	No Charge
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	96.00
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	96.00
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount. Note: Customization may be required if ordered along with an Omni.</i>	529.00



## Accessories

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107189-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Mount on left armrest: _____ Mount on right armrest: _____</i>	218.00
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	104.00
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	294.00
<input type="checkbox"/> I105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	203.00
<input type="checkbox"/> I105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	203.00
<input type="checkbox"/> I105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	203.00
<input type="checkbox"/> I105262-99-0	Joystick Cutout Tray - 16" W <i>Must also select PN I10556.</i>	492.00
<input type="checkbox"/> I105263-99-0	Joystick Cutout Tray - 18" W <i>Must also select PN I10556.</i>	492.00
<input type="checkbox"/> I10843	Small Angle Adjustable Tray <i>Must also select PN I10556. Mount on Left: _____ Mount on Right: _____</i>	492.00
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	144.00
<input type="checkbox"/> I108230-99-0	Permobil Adjustable Tray - 13"W	362.00
<input type="checkbox"/> I108228-99-0	Permobil Adjustable Tray - 15"W	362.00
<input type="checkbox"/> I108227-99-0	Permobil Adjustable Tray - 17"W	362.00
<input type="checkbox"/> I108229-99-0	Permobil Adjustable Tray - 19"W	362.00
<input type="checkbox"/> I108232-99-0	Permobil Tray Hardware Right Side Adjustable	263.00
<input type="checkbox"/> I108233-99-0	Permobil Tray Hardware Left Side Adjustable	263.00
<input type="checkbox"/> I105224-99-0	Oxygen holder	358.00
<input type="checkbox"/> I105223-99-0	Crutch holder	372.00

## Alternative Drive Systems

## ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	3,278.00
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,599.00
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,830.00

## ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA
-------------	-------------	-----------

## Alternative Drive Systems

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA
	NOTE: Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. The Swing Away Chin Boom Assemblies are Compatible with Permobil Headrests I107470-99-0, I107469-99-0 or I107468-99-0 and must be selected with I108721-99-0.	
<input type="checkbox"/> I107946-99-0	Swing Away Chin Boom Assm for Compact Lite - Left	593.00
<input type="checkbox"/> I107948-99-0	Swing Away Chin Boom Assm for Compact Lite - Right	593.00
<input type="checkbox"/> I107857-99-0	Swing Away Chin Boom Assm for Mini Joystick - Left	593.00
<input type="checkbox"/> I107947-99-0	Swing Away Chin Boom Assm for Mini Joystick - Right	593.00
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest Mount in front of Left Arm: _____ Mount in front of Right Arm: _____	311.00
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____	629.00
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.	1,007.00
<input type="checkbox"/> I10752	Haims Harness for Chin Control Mount, Large For neck width up to 19cm.	719.00
<input type="checkbox"/> I10753	Haims Harness for Chin Control Mount, Small For neck width up to 14cm.	719.00

### ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode & On/Off function. Switch operation force = 130 g Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____	298.00
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode & On/Off function. Switch operation force = 130 g	543.00

### SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10573	Sip & Puff System for Omni Requires Omni Control Module. Includes breath tube kit.	1,923.00
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop Requires Omni Control Module. Includes breath tube kit.	2,559.00
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.	1,062.00
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff Required when selecting PN I10573, I10574 or I106807-99-0.	311.00

### PROPORTIONAL FOOT CONTROL:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10767	Proportional Foot Control Not available with Active Reach Packages. Mount on Left: _____ Mount on Right: _____	3,376.00
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control Required when selecting PN I10767.	311.00

## Total Control Head Array

## Total Control Head Array

### TOTAL CONTROL HEAD ARRAY:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Must select hardware below to complete assembly. ESP recommended. Choose Omni Control Module separately.</i>	5,933.00
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes &amp; Comments section. Choose Omni Control Module separately.</i>	No Charge
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):		
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links	330.00
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	545.00
SELECT LEFT ARM PACKAGE (OPTIONAL)		
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor	449.00
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch	299.00
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch)	183.00
SELECT RIGHT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor	449.00
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch	299.00
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch)	183.00

### ASL FUSION HEAD ARRAY:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I109582-99-0	ASL Fusion Head Array (available from 4/1/2022) <i>Proportional and digital head array includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	8,515.00
<input type="checkbox"/> I11580	Swing Away Hardware for Lateral Sensor - Left	425.00
<input type="checkbox"/> I11581	Swing Away Hardware for Lateral Sensor - Right	425.00
<input type="checkbox"/> I11582	Adj/Removable Headrest Mounting Hardware	425.00
<input type="checkbox"/> I109583-99-0	ASL Fusion Display and Programmer with Mount <i>Includes ASL 165 and mounting bracket.</i>	905.00

## Alt Drive Systems Setup

### CHOOSE OMNI (REQUIRED):

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108500-99-0	Omni2 Control Module for R-net - Corpus <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	2,095.00
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	305.00

### CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge

## Alt Drive Accessories

### PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	334.00
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	587.00

## Alt Drive Accessories

## PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	846.00
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	437.00
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	291.00
<input type="checkbox"/> I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	177.00
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	153.00

## OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I101769-99-0	Microlight Switch Qty: _____	192.00
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	145.00
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	160.00
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	142.00
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	672.00
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	1,141.00
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	599.00
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	898.00
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	41.00
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	45.00
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	49.00
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	55.00
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	45.00

Notes & Additional Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgment:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.